



**20" Wide  
Stainless Steel**

- WCM2025TC27
- WCM2027TC27
- WCM2025TC50
- WCM2027TC50



Fig	Model	Description	Page
01	BR8006L-01	Brake Complete, Left *	8
	BR8006R-01	Brake Complete, Right *	8
02	BR8006X-01	Rubber Tip for Brake	9
3C	AR2605L-01-CL	Detachable Left Side Full Length Arm, Complete with skirtguard, armrest & button lock, Color * **	11
	AR2605R-01-CL	Detachable Right Side Full Length Arm, Complete with skirtguard, armrest & button lock, Color * **	11
04	AR1204X-01-CL	Full Length Armrest, Color **	12
05	AR9011X-01	Button Lock	12
06	AR3049X-01	Skirtguard 19 3/8" x 6"	12
07	WH1239X-01	12 1/2" Wheel & 7/16" Bearings	5
10	BE1018X-01	7/16" Bearing for Rear Wheel	3
11	WH1240X-01	Axle and Nut for Rear Wheel	8
14	WH3106X-01	8" Fork Assembly (fork, 8" wheel, bearings & axle)	6
15	WH3109X-01	Fork with Nut	6
16	WH3107X-01	8" Wheel with Bearings	5
17	HW3108X-01	Axle and Nut for Front Wheel	8
18	BE2019X-01	Bearing for Fork Stem	4
19	BE1012X-01	Bearing for Front Wheel	4
20	UP1315X-01-CL	Back Upholstery w/Hardware **	17
21	UP0710-01-CL	Seat Upholstery w/Hardware **	17
22	UP1220X-01	Seat Attaching Hardware Set	19
23	UP1222X-01	Back Attaching Hardware Set	19

Fig	Model	Description	Page
24	FR3104L-01	Footrest Complete, Left Side *	13
	FR3104R-01	Footrest Complete, Right Side *	13
25	LR3104L-01-CL	Legrest Complete Left, Color * **	14
	LR3104R-01-CL	Legrest Complete Right, Color * **	14
26	LR3108X-01-CL	Calf Pad, Color **	14
27	AC1506X-01	Heel Loop with Hardware	23
28	CA8016X-01	Bumper	9,16
29	CA8007X-01	Handle Grip	9
30	CA8006X-01	Kick Tube Cap	9
31	FR3109L-01	Foot Plate w/Heel Loop, Left *	16
	FR3109R-01	Foot Plate w/Heel Loop, Right *	16
32	FR3108X-01	Foot Plate Only	16
33	LR7012X-01	Tension "O" Ring for Foot Plate	16
34	LR7013X-01	Screw for Foot Plate	16
35	CA3026X-01	Seat Guide	10
36	CA3028X-01	Journal Cap for Fork Stem	9
38	CA8019X-01	Pivot Link and Bracket & Hdw	22
39	LR7000L-01	Lock Lever Release Kit, Left *	15
	LR7000R-01	Lock Lever Release Kit, Right *	15
40	LR7001L-01	Spring Only for Lock Lever, Left *	15
	LR7001R-01	Spring Only for Lock Lever, Right *	15
41	FM2003X-01	Cross Brace Set	21
42	LR7020L-01	Adjustment Rod Complete, Left*	15
	LR7020R-01	Adjustment Rod Complete, Right*	15
43	LR7003X-01	Calf Pad Bracket	15
44	LR7002X-01	Calf Pad Hinge Plate	15
46	UP1212X-01-CL	Handle Closing Straps **	20
50	AC1514X-01	Safety Belt, for 20" Wide Chairs	23

\* Determine Left or Right when sitting in the wheelchair  
 \*\* Upholstery Colors: 22-Burgundy, 50-Green, 02-Dk. Blue