



Stainless Steel

STM1550-01
STM1551-01
STM1552-01
STM1553-01
STM1585-01
STM1586-01
STM1587-01
STM1588-01

Fig	Model	Description
01	STM2634-01	Attaching Bolts for Head Support, Set of 2
02	STM2602-01	Red Stretcher Bumper
03	STM2603-01	Latch Strip w/ Attaching Hardware
04	STM2604-01	Back w/ Attaching Hardware
05	STM2605-01	Back Support
06	STM2633-01	Release Knob Complete for Side Rail
07	STM2455-01	Side Rail Complete, Serial Number STR-0711 & Up Lower Serial number, Please Call
08	STM2635-01	Plastic Bushing Set, for Side Rail Vertical Tube
09	STM2624-01	Frame Cross Tube Support for Fowler Crank Position Bolt
11	STM2611-01	Stretcher Main Frame, Left *
12	STM2611-02	Stretcher Main Frame, Right *
13	STM2613-01	Frame Cross Tube Support
14	STM2815-01	U Shape Frame Support Brace, with Hardware
15	STM2822-01	Railing Screws set for 1 Rail

Fig	Model	Description
16	STM2616-01	Utility Shelf
17	STM2617-01	Chain w/ Screw
18	WH3107X-01	8" Caster Wheel w/ Bearings
19	HW3108X-01	8" Caster Axle
20	BE1012X-01	8" Caster Bearings
24	STM2650-01	Adapter for Stretcher Leg, w/ Attaching Hardware
25	STM2615-01	Attaching Bolt and Nut for Caster
26	STM2651-01	5" Caster w/ Brake, Complete
27	STM2652-01	8" Caster w/ Brake and Pedal, Complete
28	STM2653-01	8" Caster Complete
30	STM1556-01	IV Pole for Stretcher
31	STM1559-01	O2 Tank Holder for E Cylinder Tanks
32	STM2460-01	Complete Crank Assembly
	STM1555-01-CL	Replacement Stretcher Pad **
	STM1558-01	Safety Straps, Set of 3
	STM1157-01-CL	Rail Pads, Pair **

* Determine Left or Right when sitting in the wheelchair

** Upholstery Colors: 22-Burgundy, 50-Green, 02-Dk. Blue