




CREATIVE LIVING MEDICAL //

MRI SAFETY SOLUTIONS AND NON-MAGNETIC EQUIPMENT

MARCH

2025



MRI DOCTOR STOOL PG.16



MRI FOLDING WALKER PG.14



MRI WHEELCHAIR PG.1



SCAN QR CODE OR VISIT WWW.CLMED.COM

2025
PRODUCT GUIDE //

www.clmed.com



ourterms//

Office Hours: Monday thru Friday 7:00 am to 3:30 pm (CST)

Terms: Our terms are NET 30 DAYS on all new and current accounts. Payment required in Advance if not current, even if the order is complete or not. Shipments will be held if an account contains invoices over 30 days past due. The account will go on permanent Payment in Advance Status when there is a history of past dues. Third party collection fees and legal fees must be paid in full before an account can be reopened.

Prices: All prices are subject to change without notice. Prices reflect single unit cost, except where noted.

Freight: Freight will be prepaid and added to your invoice, unless a freight account is provided. Freight charges are part of your invoice and are due accordingly. Freight charges are not discountable.

Proof Of Delivery: We will provide proof of delivery at your request. The fee for proof of delivery is \$15 if the customer received the package, no charge if it was lost in delivery. Be absolutely sure you did not receive the package before calling for proof of delivery.

Returns: No item will be accepted for return after 30 days from date of purchase. All returns must be labeled with a return authorization number which can be obtained from our customer service department. **Claims for shortage, breakage or any other reason must be made within 48 hours of receipt. Shipping containers must be kept if shipment is damaged or short. All returns are subject up to a 35% restocking fee. The customer is responsible for all freight charges incurred.**

Cancellation Charges: Orders may be canceled without charge providing your request is received prior to shipment. If the merchandise has been shipped there will be up to a 35 % restocking fee, plus all freight charges incurred. Please get cancellation confirmation to avoid these charges.

Other Order Information: We reserve the right to split shipments. We sell on first come first serve basis. Where available, substitutions will be shipped if the product ordered is not available or discontinued.

ourindex//

CLIPBOARDS/30.

COIL DRAPING/46-47.

DIAGNOSTICS/45.

DOCTOR STOOLS, OVERBEAD TABLES, INJECTION STANDS, AND MAYO STANDS/16-18.

EARPLUGS/32.

FIRE EXTINGUISHERS/49.

FURNITURE/20-21.

HAMPERS/24.

LADDERS/49.

MAINTENANCE/50-51.

MRI ACCESSORIES/30-31.

MRI BIOPSY/48.

MRI PROCEDURES/44-46.

MRI SHELVING/18-19.

MRI TRANSFER BOARDS/52-53

NON-FERROMAGNETIC EQUIPMENT/25-29.

NON-MAGNETIC ACCESSORIES/36-37.

OXYGEN & I.V./29, 43-45.

PATIENT COMFORT/36-37.

PPE/37.

PVC CARTS/23.

SANDBAGS/31.

SIGNS & LABELS/38-42.

STEP STOOLS/15.

STRETCHERS/12-13, 26.

SUCTION/43.

TOOLS/49.

UTILITY TABLES/CARTS/22-23.

WALKERS/14, 28.

WHEELCHAIR PARTS/ 7-11.

WHEELCHAIRS & ACCESSORIES/1-7.



18" WIDE WHEELCHAIRS

See pages 7-11 for wheelchair parts and accessories



18" WIDE WHEELCHAIR WITH FIXED FOOTRESTS

Fixed Footrests

BURGUNDY MODEL # WCM1825FF27

GREEN MODEL # WCM1825FF50

DARK BLUE MODEL # WCM1825FF02

with Removable Full Length Padded Armrests

SPECIFICATIONS:

Seat Width:	18"
Overall Width:	26"
Overall Depth:	44"
w/ footrests	
Seat Depth:	16"
Overall Height:	35"
Back Height:	18"
Front Casters:	8" x 1" solid
Rear Wheel:	24" solid
Wheel Locks:	push to lock
Weight:	44 lbs.

250 LB WEIGHT CAPACITY



BACK VIEW

Ship Weight: 54 lbs
Ship Dimensions: 12" x 38" x 37"

WARRANTY:

5 Years on Frame & Crossbraces
1 Year on Moving Parts
90 Days on Soft Goods



18" WIDE WHEELCHAIRS

See page 7-11 for wheelchair parts and accessories

18" WIDE WHEELCHAIR WITH SWINGAWAY FOOTRESTS

Swingaway Detachable Footrests

BURGUNDY MODEL # WCM1827SD27

GREEN MODEL # WCM1827SD50

DARK BLUE MODEL # WCM1827SD02

with Flip-Up Desk Length Padded Arms

Ship Weight: 53 lbs

Ship Dimensions: 12" x 32" x 37"

SPECIFICATIONS:

Seat Width:	18"
Overall Width:	26"
Overall Depth:	44"
w/ footrests	
Seat Depth:	16"
Overall Height:	35"
Back Height:	18"
Front Casters:	8" x 1" solid
Rear Wheel:	24" solid
Wheel Locks:	push to lock
Weight:	44 lbs.

250 LB WEIGHT CAPACITY

18" WIDE WHEELCHAIR WITH SWINGAWAY LEGRESTS

Swingaway Detachable Legrests

BURGUNDY MODEL # WCM1877SD27

GREEN MODEL # WCM1877SD50

DARK BLUE MODEL # WCM1877SD02

with Flip-Up Desk Length Arms

Ship Weight: 53 lbs

Ship Dimensions: 12" x 32" x 37"



BACK VIEW

NEW - Parts breakdowns of MRI Wheelchairs are available in PDF format by e-mail, fax, mail and on www.CLMed.com



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

WARRANTY:

5 Years on Frame & Crossbraces
1 Year on Moving Parts
90 Days on Soft Goods



20" WID **TRANSPORT CHAIRS**

See page 7-11 for wheelchair parts and accessories



20" WIDE TRANSPORT CHAIR

Detachable Footrests

BURGUNDY MODEL # **WCM2025TC27**

GREEN MODEL # **WCM2025TC50**

DARK BLUE MODEL # **WCM2025TC02**

with Removable Full Length
Padded Armrests

Ship Weight: 48 lbs

Ship Dimensions: 12" x 29" x 38"

SPECIFICATIONS:

- Seat Width: 20"
- Overall Width: 25"
- Overall Depth: 27"
w/o front rigging
- Seat Depth: 16"
- Overall Height: 35.75"
- Back Height: 18"
- Front Casters: 8" x 1" solid
- Rear Wheel: 12.50" x 2.50"
- Wheel Locks: pull to lock
- Weight: 45 lbs.

w/o Front Rigging

250 LB WEIGHT CAPACITY

WARRANTY:

- 5 Years on Frame & Crossbraces
- 1 Year on Moving Parts
- 90 Days on Soft Goods



BACK VIEW



20" WIDE TRANSPORT CHAIR

Detachable Elevating Legrests

BURGUNDY MODEL # **WCM2027TC27**

GREEN MODEL # **WCM2027TC50**

DARK BLUE MODEL # **WCM2027TC02**

with Removable Full Length
Padded Armrests

Ship Weight: 48 lbs

Ship Dimensions: 12" x 33" x 37"

250 LB WEIGHT CAPACITY



20" WIDE WHEELCHAIRS

See page 7-11 for wheelchair parts and accessories

20" WIDE WHEELCHAIR

Detachable Footrests

BURGUNDY MODEL # WCM2025SD27

GREEN MODEL # WCM2025SD50

DARK BLUE MODEL # WCM2025SD02

with Removable Full Length

Padded Armrests

Ship Weight: 58 lbs

Ship Dimensions: 12" x 33" x 36"

SPECIFICATIONS:

Seat Width: 20"

Overall Width: 28"

Overall Depth: 33"

w/o front rigging

Seat Depth: 16"

Overall Height: 38"

Back Height: 18"

Front Casters: 8" x 1" solid

Rear Wheel: 24" solid

Wheel Locks: push to lock

Weight: 46 lbs.

w/o Front Rigging

250 LB WEIGHT CAPACITY

WARRANTY:

5 Years on Frame & Crossbraces

1 Year on Moving Parts

90 Days on Soft Goods



BACK VIEW



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



20" WIDE WHEELCHAIR

Detachable Elevating Legrests

BURGUNDY MODEL # WCM2075SD27

GREEN MODEL # WCM2075SD50

DARK BLUE MODEL # WCM2075SD02

with Removable Full Length

Padded Armrests

Ship Weight: 58 lbs

Ship Dimensions: 12" x 37" x 38"

250 LB WEIGHT CAPACITY

NEW - Parts breakdowns of MRI Wheelchairs are available in PDF format by e-mail, fax, mail and on www.CLMed.com

22" WIDE WHEELCHAIRS 

See page 7-11 for wheelchair parts and accessories



22" WIDE WHEELCHAIR

Detachable Footrests

BURGUNDY MODEL # WCM2226SD27

GREEN MODEL # WCM2226SD50

DARK BLUE MODEL # WCM2226SD02

with Removable Desk Length

Padded Armrests

Ship Weight: 64 lbs

Ship Dimensions: 12" x 35" x 39"

SPECIFICATIONS:

Seat Width: 22"

Overall Width: 30"

Overall Depth: 33"

w/o front rigging

Seat Depth: 18"

Overall Height: 38"

Back Height: 18"

Front Casters: 8" x 1" solid

Rear Wheel: 24" solid

Wheel Locks: push to lock

Weight: 46 lbs.

w/o Front Rigging

350 LB WEIGHT CAPACITY



BACK VIEW



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



22" WIDE WHEELCHAIR

Detachable Elevating Legrests

BURGUNDY MODEL # WCM2276SD27

GREEN MODEL # WCM2276SD50

DARK BLUE MODEL # WCM2276SD02

with Removable Desk Length

Padded Armrests

Ship Weight: 64 lbs

Ship Dimensions: 13" x 40" x 35"

350 LB WEIGHT CAPACITY

NEW - Parts breakdowns of MRI Wheelchairs are available in PDF format by e-mail, fax, mail and on www.CLMed.com

24" WIDE WHEELCHAIRS

See page 7-11 for wheelchair parts and accessories

24" WIDE WHEELCHAIR

Detachable Footrests

BURGUNDY MODEL # WCM2424SD27

GREEN MODEL # WCM2424SD50

DARK BLUE MODEL # WCM2424SD02

with Removable Desk Length

Padded Armrests

Ship Weight: 65 lbs

Ship Dimensions: 12" x 39" x 34"

SPECIFICATIONS:

Seat Width: 24"

Overall Width: 32"

Overall Depth: 33"

w/o front rigging

Seat Depth: 18"

Overall Height: 38"

Back Height: 18"

Front Casters: 8" x 1" solid

Rear Wheel: 24" solid

Wheel Locks: push to lock

Weight: 61 lbs.

w/o Front Rigging

350 LB WEIGHT CAPACITY

WARRANTY:

5 Years on Frame & Crossbraces

1 Year on Moving Parts

90 Days on Soft Goods



BACK VIEW



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



24" WIDE WHEELCHAIR

Detachable Elevating Legrests

BURGUNDY MODEL # WCM2474SD27

GREEN MODEL # WCM2474SD50

DARK BLUE MODEL # WCM2474SD02

with Removable Desk Length

Padded Armrests

Ship Weight: 65 lbs

Ship Dimensions: 13" x 40" x 35"

350 LB WEIGHT CAPACITY

NEW - Parts breakdowns of MRI Wheelchairs are available in PDF format by e-mail, fax, mail and on www.CLMed.com



24" WIDE BARIATRIC HD CHAIRS

See page 7-11 for wheelchair parts and accessories



24" WIDE WHEELCHAIR

Detachable Footrests
BURGUNDY MODEL # WCM2425SD27
GREEN MODEL # WCM2425SD50
DARK BLUE MODEL # WCM2425SD02

NOTE: (OVERALL WIDTH IS 33")

with 20" Seat Depth,
 Removable Desk Length
 Padded Armrests
 Solid Seat Insert

24" WIDE WHEELCHAIR

Detachable Elevating Legrests
BURGUNDY MODEL # WCM2475SD27
GREEN MODEL # WCM2475SD50
DARK BLUE MODEL # WCM2475SD02

NOTE: (OVERALL WIDTH IS 33")

with 20" Seat Depth,
 Removable Desk Length
 Padded Armrests
 Solid Seat Insert

Ship Weight: 66 lbs
Ship Dimensions: 13" x 37" x 40"

SPECIFICATIONS:

Seat Width:	24"
Overall Width:	33"
Seat Width:	24"
Overall Depth:	36"
w/o front rigging	
Seat Depth:	20"
Overall Height:	38"
Back Height:	18"
Front Casters:	8" x 2" solid
Rear Wheel:	24" Pneumatic
Wheel Locks:	push to lock
Weight:	58 lbs.
w/o Front Rigging	

550 LB WEIGHT CAPACITY

WARRANTY:

5 Years on Frame & Crossbraces
 1 Year on Moving Parts
 90 Days on Soft Goods



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

NEW - Parts breakdowns of MRI Wheelchairs are available in PDF format by e-mail, fax, mail and on www.CLMed.com

26" WIDE BARIATRIC HD CHAIRS



See page 7-11 for wheelchair parts and accessories

26" WIDE WHEELCHAIR

Detachable Footrests
BURGUNDY MODEL # WCM2626SD27
GREEN MODEL # WCM2626SD50
DARK BLUE MODEL # WCM2626SD02

NOTE: (OVERALL WIDTH IS 35")

with 20" Seat Depth,
 Removable Desk Length
 Padded Armrests
 Solid Seat Insert

26" WIDE WHEELCHAIR

Detachable Elevating Legrests
BURGUNDY MODEL # WCM2676SD27
GREEN MODEL # WCM2676SD50
DARK BLUE MODEL # WCM2676SD02

NOTE: (OVERALL WIDTH IS 35")

with 20" Seat Depth,
 Removable Desk Length
 Padded Armrests
 Solid Seat Insert

Ship Weight: 66 lbs
Ship Dimensions: 14" x 37" x 40"

SPECIFICATIONS:

Seat Width:	26"
Overall Width:	35"
Overall Depth:	36"
w/o front rigging	
Seat Depth:	20"
Overall Height:	38"
Back Height:	18"
Front Casters:	8" x 2" solid
Rear Wheel:	24" Pneumatic
Wheel Locks:	push to lock
Weight:	61 lbs.

550 LB WEIGHT CAPACITY



BACK VIEW



WARRANTY:

5 Years on Frame & Crossbraces
 1 Year on Moving Parts
 90 Days on Soft Goods

WHEELCHAIR ACCESSORIES 

STAINLESS STEEL ANTI THEFT BARS

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
AC1520X-01	for 18" Wide Wheelchair	24 lbs	20" x 5" x 46"
AC1521X-01	for 20" Wide Wheelchair	24 lbs	20" x 5" x 46"
AC1522X-01	for 22" Wide Wheelchair	24 lbs	20" x 5" x 46"
AC1523X-01	for 24" Wide Wheelchair	40 lbs	24" x 7" x 46"



ANTI-TIPPERS

MODEL # **AC1507X-01**

Ship Weight: 3 lbs

Ship Dimensions: 5" x 11" x 11"



I.V. POLE

MODEL # **ACM1500-01**

Ship Weight: 12 lbs

Ship Dimensions: 8" x 8" x 36"



I.V. & OXYGEN TANK HOLDER

MODEL # **ACM1501-01**

Ship Weight: 12 lbs

Ship Dimensions: 8" x 8" x 36"



NON-MAGNETIC SAFETY PADLOCK

Plastic shackle helps prevent electrocution and arc flashing, also prevents electrical current from shackle to key
All components are non-sparking and non-magnetic - plastic, copper, brass
1½" wide 1¾" tall Xenoy green or red body, plastic shackle w/1½" clearance

Designed exclusively for Lockout/Tagout applications
Durable, lightweight, non-conductive Xenoy lock body

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
PLM1000-01	Green	2 Lbs	7" x 4" x 4"
PLM1001-01	Red	2 Lbs	7" x 4" x 4"

ANTI-THEFT CLAMP

Anti-Theft Clamp, Set of 2

MODEL # **CA8010X-01**

Ship Weight: 5 lbs

Ship Dimensions: 7" x 4" x 4"



OXYGEN TANK HOLDER

MODEL # **ACM1502-01**

Ship Weight: 10 lbs

Ship Dimensions: 10" x 10" x 19"



STENCIL COMPANY NAME OR OTHER INFORMATION

MODEL # **ACM1503-01**

Ship Weight: 5 lbs

Ship Dimensions: 4" x 12" x 19"



Option available for no weight capacity stenciling.

MRI GUIDES, GRIPS, & CAPS

RAIL GLIDE

MODEL # **CA3029X-01**
White Plastic for MRI Wheelchairs



SEAT GUIDE

MODEL # **CA3026X-01**
Black Plastic for MRI Wheelchairs



SEAT GUIDE WITH ARM SOCKET

MODEL # **CA3027X-01**
Black Plastic for MRI Wheelchairs



BRAKE/WHEEL LOCK GRIP

MODEL # **BR8007X-01**
Black Rubber for MRI Wheelchairs
Grip Slot Size: 5/8" x 1/8"



HAND GRIP

MODEL # **CA8007X-01**
Black for MRI Wheelchairs



KICK TUBE BACK

MODEL # **CA8006X-01**
Black Rubber for MRI Wheelchairs
Size: Fits 7/8" Tubing



JOURNAL CAP

MODEL # **CA3028X-01**
Black Plastic for MRI Wheelchairs
Cap Size: 1 13/64" I.D.



HUB CAP

MODEL # **CA3024X-01**
Black Plastic for 18" to 24" Wide
Standard MRI Wheelchairs



HUB CAP

MODEL # **CA3025X-01**
Black Plastic for 24" to 26"
Heavy Duty MRI Wheelchairs



RUBBER BUMPER

MODEL # **CA8016X-01**
for Footrest/Legrest of MRI Wheelchairs
Fits: 3/4" I.D. Tubing



MRI ARMREST ASSEMBLY

COMPLETE DESK LENGTH FLIP BACK ARM ASSEMBLY

MODEL # **AR2604SD-01-CL**
With Padded Armrest, for 18" Wide
Wheelchairs, Choose Side, Color



COMPLETE FULL LENGTH DETACHABLE ARM ASSEMBLY

MODEL # **AR2605SD-01-CL**
With Padded Armrest, for 18" & 20" Wide
Wheelchairs, Choose Side, Color



COMPLETE DESK LENGTH DETACHABLE ARM ASSEMBLY

MODEL # **AR2607SD-01-CL**
With Padded Armrest, for 22" & 24" Wide
Wheelchairs, Choose Side, Color



COMPLETE FULL LENGTH DETACHABLE ARM ASSEMBLY

MODEL # **AR2606SD-01-CL**
With Padded Armrest, for 24" &
26" Wide Heavy Duty Wheelchairs,
Choose Side, Color



RELEASE LEVER ASSEMBLY

MODEL # **CA8008X-01**
Black Plastic for Flip Back Arms on MRI Wheelchairs



BUTTON LOCK

MODEL # **AR9011X-01**
Stainless Steel, for Removable
Arms on MRI Wheelchairs
Button Size: 3/8" Long



DESK LENGTH PADDED ARMREST

MODEL # **AR2204X-01-CL**
Hole Spacing: 3 5/8" with Brass Inserts
Overall Length: 10 1/4" Choose Color



FULL LENGTH PADDED ARMREST

MODEL # **AR1204X-01-CL**
Hole Spacing: 9" with Brass Inserts
Overall Length: 14 1/2" Choose Color



SKIRT GUARD FOR DESK LENGTH REMOVABLE ARMS

MODEL # **AR3048X-01**
Size: 13 1/4" Long x 6" High
Stainless Steel



SKIRT GUARD FOR FULL LENGTH REMOVABLE ARMS

MODEL # **AR3049X-01**
Size: 19 1/4" Long x 6" High
Stainless Steel



MRI UPHOLSTERY

CHOOSE COLOR: **BURGUNDY** - 22 | **GREEN** - 50 | **DARK BLUE** - 02



SEAT 18" WIDE X 16" DEEP
MODEL # **UP0310X-01-CL**
Hole Spacing (from front to back): 2" to 4" to 4" to 4"



SEAT 20" WIDE X 16" DEEP
MODEL # **UP0710X-01-CL**
Hole Spacing (from front to back): 2" to 4" to 4" to 4"



SEAT 22" WIDE X 18" DEEP
Model # **UP0410X-01-CL**
Hole Spacing (from front to back): 3" to 4" to 4" to 4"



SEAT 24" WIDE X 18" DEEP
MODEL # **UP0610X-01-CL**
Hole Spacing (from front to back): 3" to 4" to 4" to 4"



SOLID SEAT
24" WIDE X 20" DEEP
MODEL # **UP0810X-01-CL**
1 1/4" foam, 1/2" ply-wood, includes metal straps. For 24" HD Wheelchairs



SOLID SEAT
26" WIDE X 20" DEEP
MODEL # **UP0510X-01-CL**
1 1/4" foam, 1/2" ply-wood, includes metal straps. For 26" HD Wheelchairs



BACK 18" WIDE X 18" TALL
MODEL # **UP1310X-01-CL**
Hole Spacing (from top to bottom): 3/4" to 5 1/2" to 5 1/2" to 5 1/2"



BACK 20" WIDE X 18" TALL
MODEL # **UP1315X-01-CL**
Hole Spacing (from top to bottom): 3" to 5 1/2" to 5 1/2" to 5 1/2"



BACK 22" WIDE X 18" TALL
MODEL # **UP1312X-01-CL**
Hole Spacing (from top to bottom): 3" to 5 1/2" to 5 1/2" to 5 1/2"



BACK 24" WIDE X 19" TALL
MODEL # **UP1314X-01-CL**
Hole Spacing (from top to bottom): 3" to 5 1/2" to 5 1/2" to 5 1/2"



BACK 24" WIDE X 19" TALL
MODEL # **UP1318X-01-CL**
Slip-On Back
Hole Spacing (from top to bottom): 7/8"



BACK 26" WIDE X 19" TALL
MODEL # **UP1313X-01-CL**
Slip-On Back
Hole Spacing (from top to bottom): 7/8"

MRI WHEELCHAIR BEARINGS



1/2" ID X 1-1/16" OD
MODEL # **BE2019X-01**
with Flange for Fork / Caster Stem, for MRI Compatible Wheelchairs



5/16" ID X 29/32" OD
MODEL # **BE1012X-01**
for Front 8" x 1" Wheel, for 18" to 24" Wide Standard MRI Wheelchairs



5/16" ID X 7/8" OD
MODEL # **BE1017X-01**
for Front 8" x 2" Wheel, for 24" and 26" Wide Heavy Duty Wheelchairs



7/16" ID X 1 1/8" OD
MODEL # **BE1014X-01**
with 1 3/8" Flange, for 24" Wheel of MRI compatible 18" and 20" Wide Wheelchairs



5/8" ID X 1 1/4" OD
MODEL # **BE1013X-01**
with 1 3/8" Flange, for 24" Wheel of MRI compatible 22" and 24" Wide Wheelchairs



5/8" ID X 1 1/4" OD
MODEL # **BE1015X-01**
with 1 3/8" Double Flange, for 24" Wheel of MRI compatible 24" and 26" Wide Heavy Duty Wheelchairs



5/8" ID X 1 1/4" OD
MODEL # **BE1016X-01**
Ceramic Ball Bearing, with 1 3/8" Flange, for 24" Wheel of MRI compatible Heavy Duty Wheelchairs

MRI FRONT RIGGING

CLAMP-ON STYLE MRI REPLACEMENT EXTENSION TUBE ASSEMBLY & FOOTREST

MODEL # **FR8808L-01** (Left)
MODEL # **FR8808R-01** (Right)
for 18" Wide Fixed Chair, Choose Side



DETACHABLE SWING-AWAY FOOTREST WITH HEEL LOOPS

MODEL # **FR3104L-01** (Left)
MODEL # **FR3104R-01** (Right)
for 18" & 20" Wide Wheelchairs Choose Side
MODEL # **FR3105L-01**
MODEL # **FR3105R-01**
for 22" & 24" Wide Wheelchairs Choose Side
MODEL # **FR3107L-01**
MODEL # **FR3107R-01**
for 24" & 26" Wide Heavy Duty Wheelchairs, Choose Side



ALUMINUM DETACHABLE SWING-AWAY FOOTREST WITH HEEL LOOPS

MODEL # **FR3201L-05** (Left)
MODEL # **FR3201R-05** (Right)
for 20" & 22" Wheelchairs
MODEL # **FR3202L-05** (Left)
MODEL # **FR3202R-05** (Right)
for 24" Wheelchairs



DETACHABLE SWING-AWAY LEGREST WITH CALF PAD

MODEL # **LR3104L-01-CL**
MODEL # **LR3104R-01-CL**
for 18" to 24" Wide Standard Wheelchairs,
Choose Side, Choose Color
MODEL # **FR3107L-01-CL**
MODEL # **FR3107R-01-CL**
for 24" to 26" Wide Heavy Duty
Wheelchairs, Choose Side, Choose Color



PADDED CALF PAD

Hole Spacing: 2", Universal, Choose Color
MODEL # **LR3108X-01-CL**
Burgundy
Green
Dark Blue



FOOTPLATE

MODEL # **FR3108X-01**
for 18" to 24" Wide Standard
Wheelchairs, Silver Aluminum
Heel loops not included
MODEL # **FR3111X-01**
for 24" to 26" Wide Heavy Duty
Wheelchairs, Silver Aluminum
Heel loops not included



HEEL LOOP FOR FOOTREST

Hardware Included
MODEL # **AC1506X-01**
for 18" & 24" Wide Wheelchairs
Fits Left or Right



MODEL # **AC1518L-01**
MODEL # **AC1518R-01**
for 24" & 26" Wide Heavy Duty
Wheelchairs, Choose Side



SILVER ALUMINUM FOOTPLATE WITH HEEL LOOP

MODEL # **FR3109L-01**
MODEL # **FR3109R-01**
for 18" to 24" Wide Wheelchairs
Choose Side



MODEL # **FR3112L-01**
MODEL # **FR3112R-01**
for 24" & 26" Wide Heavy Duty
Wheelchairs, Choose Side



LOCKING KIT FOR FOOTREST OR LEGREST

MODEL # **CA8011X-01**
Black Plastic, Sold in Pairs
Size: Fits 1/4" Front Rigging Pin
Allen Wrench Included



FOOTREST TENSION "O" RING

MODEL # **LR7012X-01**
Black Rubber, Fits all MRI Wheelchairs



BUMPER SCREW

MODEL # **LR7013X-01**
Holds bumper in and prevents
footplate from sliding off.
Fits all MRI Wheelchairs



FOOTREST / LEGREST BUMPER

MODEL # **CA8016X-01**
Holds bumper in and prevents
footplate from sliding off.
Fits all MRI Wheelchairs



FOOTREST / LEGREST ROLLER BUMPER

MODEL # **CA8004X-01**
for Heavy Duty Front Rigging ONLY
Inserts in the end of the footplate tube
Black Rubber and Plastic





MRI ALUMINUM BRAKE
MODEL # **BR8008L-01**
MODEL # **BR8008R-01**
for Detachable Arm Wheelchair
Tube Mount, Choose Side

MRI BRAKES / WHEEL LOCKS



BRAKE GRIPS
Model # **BR8007X-01**
Black Rubber, for MRI Wheelchairs
Grip Slot Size: 5/8" x 1/8"



MRI STAINLESS BRAKE
MODEL # **BR8006L-01**
MODEL # **BR8006R-01**
for Detachable Arm Wheelchair
Tube Mount, Choose Side

MRI WHEELS



FRONT FORK WITH NUT
MODEL # **WH9109X-01**
for Standard and Heavy
Duty Wheelchairs



8" X 1" FRONT WHEEL
MODEL # **WH3107X-01**
for 18" to 24" Wide
Standard Wheelchairs



8" X 2" FRONT WHEEL
MODEL # **WH3109X-01**
for 24" and 26" Wide Heavy
Duty Wheelchairs



8" X 1" WHEEL ASSEMBLY
MODEL # **WH3106X-01**
for 18" to 24" Wide
Standard Wheelchairs



8" X 2" WHEEL ASSEMBLY
MODEL # **WH3108X-01**
for 24" to 26" Wide Heavy
Duty Wheelchairs



BLACK PLASTIC HUB CAP
MODEL # **CA3024X-01**
18" to 24" Wide Standard
Wheelchairs



BLACK PLASTIC HUB CAP
MODEL # **CA3025X-01**
24" to 26" Wide Heavy
Duty Wheelchairs



AXLE AND NUT
MODEL # **HW6108X-01**
for 8" Front Wheel on MRI Wheelchairs
Axle Size: 5/16" Diameter x 3 1/4" Long



REAR AXLE AND NUT
MODEL # **WH6109X-01**
for 18" & 20" Stainless Steel Wheelchairs
Axle Size: 7/16" Diameter x 5" Long



REAR AXLE AND NUT
MODEL # **WH6110X-01**
for 22", 24" and 26" Standard and
HD Stainless Steel Wheelchairs
Axle Size: 5/8" Diameter x 5" Long



24" REAR WHEELS COMPLETE
MODEL # **WH1224X-01**
18" to 20" Wide, with 7/16" Bearings
MODEL # **WH1225X-01**
22" to 24" Wide, with 5/8" Bearings



HAND RIM SCREWS AND SPACERS
MODEL # **HW9213X-01**
for MRI Wheelchairs, includes 8
screws and 8 plastic spacers



**VICE TIRE SETTER FOR
SNAP-ON TIRES**
MODEL # **CA8012X-01**
Wheel Sizes: 8", 20", 22" 24"
THIS ITEM IS MAGNETIC

MRI FIXED **STRETCHER** 

- Stainless steel construction
- Drop down side rails
- Grounding chain
- Storage tray
- Full surround rubber bumper
- Adjustable headrest
- 3" Pad

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
STM1550-01	Fixed Gurney, 5" Total Lock Casters, Black Pad	225 lbs	48" x 40" x 87"
STM1551-01	Fixed Gurney, 2-8" Swivel & 2-8" Locking Casters, Black Pad	225 lbs	48" x 40" x 87"
STM1552-01	Fixed Gurney, 5" Total Lock Casters, Green Pad	225 lbs	48" x 40" x 87"
STM1553-01	Fixed Gurney, 2-8" Swivel & 2-8" Locking Casters, Green Pad	225 lbs	48" x 40" x 87"
STM1585-01	Fixed Gurney, 5" Casters, Black Pad, Includes Crank Assembly	225 lbs	48" x 40" x 87"
STM1586-01	Fixed Gurney, 8" Casters, Black Pad, Includes Crank Assembly	225 lbs	48" x 40" x 87"
STM1587-01	Fixed Gurney, 5" Casters, Green Pad, Includes Crank Assembly	225 lbs	48" x 40" x 87"
STM1588-01	Fixed Gurney, 8" Casters, Green Pad, Includes Crank Assembly	225 lbs	48" x 40" x 87"
STM1590-01	Fixed Gurney, 5" Total Lock Casters, Dark Blue	225 lbs	48" x 40" x 87"
STM1591-01	Fixed Gurney, 2-8" Swivel & 2-8" Swivel Locking Casters, Dark Blue	225 lbs	48" x 40" x 87"
STM1592-01	Fixed Gurney, 5" Total Lock Casters, Burgundy	225 lbs	48" x 40" x 87"
STM-1593-01	Fixed Gurney, 2-8" Swivel & 2-8" Swivel Locking Casters, Burgundy	225 lbs	48" x 40" x 87"

SPECIFICATIONS:

- Height: 30 1/4"
- Height with Rails up: 42"
- Height of Rails: 9"
- Length: 77"
- Rail Length: 49"
- Width: 30"
- Weight: 100 lbs.
- Back Adj. Range: 0-75 Degrees
- Number of Positions: 11
- Storage Tray: 1 1/2"D x 26"L x 20 1/4"W


New and Improved - now even stronger

WEIGHT CAPACITY: 600 LBS

- Ship Weight: 225 lbs
- Ship Dimensions: 40" x 48" x 84"

WARRANTY:

- 5 Years on Frame, 1 Year on Moving Parts
- 90 Days on Soft Goods



Specify Color



STM1585-01
GURNEY W/HAND CRANK



8" CASTER



STM1551-01
GURNEY W/8" CASTERS

STRETCHER **ACCESSORIES** 

REPLACEMENT PAD

- MODEL # STM1555-CL
- Fits Gurney Above
- Height: 3"
- Width: 22 1/2"
- Length: 77"
- Ship Weight: 30 lbs
- Ship Dimensions: 6" x 24" x 40"

SAFETY STRAPS

- MODEL # STM1558-01
- Fits Gurney Above
- Set of 3
- Ship Weight: 3 lbs
- Ship Dimensions: 5" x 11" x 11"



RAIL PADS

- MODEL # STM1157-CL
- Fits Gurney Above
- Ship Weight: 25 lbs
- Ship Dimensions: 14" x 14" x 24"

I.V. POLE

- MODEL # STM1556-01
- Fits Gurney Above
- Height: 32"
- Ship Weight: 12 lbs
- Ship Dimensions: 8" 8" x 36"



OXYGEN TANK HOLDER

- MODEL # STM1559-01
- Fits Gurney Above
- Ship Weight: 12 lbs
- Ship Dimensions: 8" x 8" x 36"

FERNO MRI **STRETCHER**

Model 30NM The "Standard" and the lightest lift-in cot in the industry. Ferno's Model 30 Multi-Level, Lift-in Cot served as the industry standard for many years and is still used by many ambulance companies and medical facilities. The Ferno Model 30 cot adjusts to six height levels, enabling patient transfer from most nursing home beds or hospital treatment tables. To operate the height adjustment mechanism, the EMTs must assume the patients weight prior to raising or lowering the cot-bed frame. Because the Model 30 must be lifted to be loaded into the vehicle, it can be used with ambulances that have various floor heights. The Model 30 cot will work well in station wagons, vans or modular ambulances. The telescoping head-end main frame section allows the Model 30 cot to be shortened for maneuverability around tight corners and spaces. The Model 30NM is the Non-Magnetic version for medical treatment areas where the cot must be non-magnetized.

SPECIFICATIONS:

- Max Height: 32"
- Min Height: 12"
- Max Length: 78"
- Width: 23"
- Weight: 56 lbs
- Number of Positions: 6
- Casters: All Locking
- WEIGHT CAPACITY: 400 LBS**

WARRANTY:

- 2 Years on Frame
- 90 Days on Pad & Straps



STM1454-01

CONTOUR STYLED PAD



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
STM1454-01	Gurney I.V. Pole	15 lbs	6" x 6" x 36"
STM1450-01	With Flat Style Pad	90 lbs	40" x 48" x 85"
STM1451-01	With Contoured Pad	90 lbs	40" x 48" x 85"

OTHER REPLACEMENT **STRETCHER PADS**

Replacement stretcher pads covered with Sure-Check Fusion III material. Available with square corners or angle cut corners. **Stretcher pads are Non-Returnable.** Custom sizes available



SQUARE CORNER PAD



CUT CORNERS PAD

MODEL #	DIMENSIONS	SHIP WEIGHT	SHIP L X W X H
STM4000-01	24"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4001-01	24"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4002-01	24"W x 75"L x 4"T	30 lbs	25" x 25" x 36"
STM4003-01	26"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4004-01	26"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4005-01	26"W x 75"L x 4"T	30 lbs	25" x 25" x 36"
STM4006-01	28"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4007-01	28"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4008-01	28"W x 75"L x 4"T	30 lbs	25" x 25" x 36"

MODEL #	DIMENSIONS	SHIP WEIGHT	SHIP L X W X H
STM4009-01	24"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4010-01	24"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4011-01	24"W x 75"L x 4"T	30 lbs	25" x 25" x 36"
STM4012-01	26"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4013-01	26"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4014-01	26"W x 75"L x 4"T	30 lbs	25" x 25" x 36"
STM4015-01	28"W x 75"L x 2"T	30 lbs	25" x 25" x 36"
STM4016-01	28"W x 75"L x 3"T	30 lbs	25" x 25" x 36"
STM4017-01	28"W x 75"L x 4"T	30 lbs	25" x 25" x 36"



MRI NON MAGNETIC CANES

Wooden canes are made from the finest hardwoods. The cane is kiln cured for maximum strength. Rubber tips provide additional safety.

Specifications:

1" x 36"

Available in Oak, Mahogany or Walnut

250 LB CAPACITY

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
CNM1000-01-82	Oak	8 lbs	36" x 6" x 6"
CNM1000-01-83	Mahogany	8 lbs	36" x 6" x 6"
CNM1000-01-86	Walnut	8 lbs	36" x 6" x 6"


TRANSFERBELT

MODEL # **GBM1500-01**

2" x 48"

Ship Weight: 1lbs

Ship Dimensions: 4" x 4" x 10"



MRI FOLDING WALKERS

18" Wide, 17" Deep

250 LB CAPACITY

JUNIOR

MODEL # **WAM1001-01**

Adj. Height 28" to 32"

Ship Weight: 23 lbs

Ship Dimensions: 6" x 24" x 30"

ADULT

MODEL # **WAM1000-01**

Adj. Height 32" to 36"

Ship Weight: 24 lbs

Ship Dimensions: 6" x 24" x 32"

HEAVY DUTY - EXTRA WIDE

MODEL # **WAM1002-01**

500 LB CAPACITY

24" Wide, 21" Deep

Adj. Height 32" to 38 1/2"

Ship Weight: 30 lbs

Ship Dimensions: 6" x 28" x 34"

WARRANTY

3 year Warranty

Excluding Tips on Walkers

5" ADJUSTABLE WALKER WHEELS

MODEL # **WAM1010-01**

350 LB CAPACITY

Set of 2

Ship Weight: 6 lbs

Ship Dimensions: 5" x 11" x 11"



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

MRI STEP STOOLS

Heavy Duty Stainless Steel
550 LB CAPACITY

MODEL # FSM1724-01

Step Stool without Handrail
16" Wide x 12" Depth
Step Height: 8 1/2"
Ship Weight: 18 lbs
Ship Dimensions: 14" x 14" x 18"



MODEL # FSM1725-01

Step Stool with 41" Handrail
16" Wide x 12" Depth
Step Height: 8 1/2"
Ship Weight: 48 lbs
Ship Dimensions: 15" x 15" x 42"



MRI NARROW STEP STOOLS

Heavy Duty Stainless Steel
550 LB CAPACITY

MODEL # FSM1730-01

Step Stool without Handrail
16" Wide x 26" Depth
Step Height: 8 1/2"
Weight: 20 lbs




MODEL # FSM1731-01

Step Stool with Handrail
16" L X 26" Depth
Step Height: 8 1/2"
Handrail: 47" Tall x 26" Deep
Weight: 20 lbs



MRI BARIATRIC STEP STOOLS

Heavy Duty Stainless Steel
550 LB CAPACITY



MODEL # FSM1733-01

Step Stool no Handrails
30" Wide x 12" Depth
Step Height: 8 1/2"
Ship Weight: 24 lbs
Ship Dimensions: 15" x 15" x 42"



MODEL # FSM1734-01

Step Stool with One 41" Handrail
30" Wide x 12" Depth
Step Height: 8 1/2"
Ship Weight: 36 lbs
Ship Dimensions: 15" x 15" x 42"



MODEL # FSM1729-01

Step Stool with Two 41" Handrails
30" Wide x 12" Depth
Step Height: 8 1/2"
Ship Weight: 48 lbs
Ship Dimensions: 15" x 15" x 42"

MRI DOUBLE STEP STOOL

Heavy Duty Stainless Steel
550 LB CAPACITY

MODEL # FSM1726-01

with 41" Rail.
36" Wide x 24" Depth
Step Height: 8 1/2"
Weight: 47 lbs
Assembled Ship Weight: 214 lbs
Assembled Ship Dimensions: 24" x 48" x 36"
Unassembled Ship Weight: 26 lbs
Unassembled Ship Dimensions: 4" x 25" x 50"



WARRANTY:

5 Years on Frame, 90 Days on Soft Goods



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.





CHM8038-01
w/ Casters, Back & Armrests

WARRANTY:
5 Years on Frame
1 Year on Moving Parts
90 Days on Soft Goods

BURGUNDY DOCTOR STOOLS

- 4" Padded Seat
- Stainless Steel Construction
- Five Leg Design for Maximum Stability
- Padded Armrests (on selected models)
- 2" Dual Wheel Casters (on selected models)

WEIGHT CAPACITY: 250 LBS



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

MODEL #	ADJ. HEIGHT	FOOT STYLE	BACKREST	ARMRESTS	SHIP WEIGHT	SHIP L X W X H
CHM8030-01	15"-21"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8031-01	15"-21"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8032-01	15"-21"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8036-01	16"-22"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8037-01	16"-22"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8038-01	16"-22"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8033-01	21"-27"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8034-01	21"-27"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8035-01	21"-27"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8039-01	22"-28"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8040-01	22"-28"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8041-01	22"-28"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"



CHM8032-01 w/ Rubber Tips
Back & Armrests

WARRANTY:
5 Years on Frame
1 Year on Moving Parts
90 Days on Soft Goods

GREEN DOCTOR STOOLS

- 4" Padded Seat
- Stainless Steel Construction
- Five Leg Design for Maximum Stability
- Padded Armrests (on selected models)
- 2" Dual Wheel Casters (on selected models)

WEIGHT CAPACITY: 250 LBS



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

MODEL #	ADJ. HEIGHT	FOOT STYLE	BACKREST	ARMRESTS	SHIP WEIGHT	SHIP L X W X H
CHM8030-02	15"-21"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8031-02	15"-21"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8032-02	15"-21"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8036-02	16"-22"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8037-02	16"-22"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8038-02	16"-22"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8033-02	21"-27"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8034-02	21"-27"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8035-02	21"-27"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8039-02	22"-28"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8040-02	22"-28"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8041-02	22"-28"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"



CUSTOM STENCILING FOR DOCTOR STOOLS

MODEL # CHM8914-01

Customize your doctor stool. Add your facility name, department name, or any information you would like.

DARK BLUE DOCTOR STOOLS



- 4" Padded Seat
- Stainless Steel Construction
- Five Leg Design for Maximum Stability
- Padded Armrests (on selected models)
- 2" Dual Wheel Casters (on selected models)



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

WEIGHT CAPACITY: 250 LBS

MODEL #	ADJ. HEIGHT	FOOT STYLE	BACKREST	ARMRESTS	SHIP WEIGHT	SHIP L X W X H
CHM8030-03	15"-21"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8031-03	15"-21"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8032-03	15"-21"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8036-01	16"-22"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8037-01	16"-22"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8038-01	16"-22"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8033-03	21"-27"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8034-03	21"-27"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8035-03	21"-27"	Rubber Tips	Yes	Yes	135 lbs	25" x 25" x 36"
CHM8039-01	22"-28"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8040-01	22"-28"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8041-01	22"-28"	Casters	Yes	Yes	135 lbs	25" x 25" x 36"



CHM8032-03w/ Rubber Tips
Back & Armrests

MRI STAINLESS DOCTOR STOOLS



- Solid non-slip Seat
- Stainless Steel Construction
- Five Leg Design for Maximum Stability
- 2" Dual Wheel Casters (on selected models)

WEIGHT CAPACITY: 250 LBS

WARRANTY:

- 5 Years on Frame
- 1 Year on Moving Parts
- 90 Days on Soft Goods



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



MODEL #	ADJ. HEIGHT	FOOT STYLE	BACKREST	ARMRESTS	SHIP WEIGHT	SHIP L X W X H
CHM8050-01	15"-21"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8051-01	15"-21"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8052-01	21"-27"	Rubber Tips	No	No	94 lbs	25" x 25" x 25"
CHM8053-01	21"-27"	Rubber Tips	Yes	No	135 lbs	25" x 25" x 36"
CHM8054-01	15"-21"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8055-01	15"-21"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8056-01	21"-27"	Casters	No	No	94 lbs	25" x 25" x 25"
CHM8057-01	21"-27"	Casters	Yes	No	135 lbs	25" x 25" x 36"
CHM8058-01	Fixed	Rubber Tips	No	No	48 lbs	20" x 20" x 20"

CREATIVE LIVING MEDICAL //

MRI SAFETY SOLUTIONS AND NON-MAGNETIC EQUIPMENT

PH: 888-622-6522 FAX: 218-824-6091
 1111 INDUSTRIAL PARK ROAD, BRAINERD MN, 56401

WATERBASIN

MODEL # SBM1350-01
 Stainless steel water basin with
 2" plastic twin wheels.
 36" Height, 14" water bowl
Ship Weight: 75 lbs
Ship Dimensions: 20" x 20" x 36"



WARRANTY:
 5 Years on Frame
 90 Days on Soft Goods

MAYOSTAND

MODEL # INM1402-01
 Stainless Steel
 Adjusts from 31" to 51"
 2" twin plastic wheels
Mayo Stand Tray size: 13" x 19"
Ship Weight: 75 lbs
Ship Dimensions: 20" x 20" x 36"



WARRANTY:
 5 Years on Frame
 90 Days on Soft Goods

OVERBEDTABLE

MODEL # INM1560-01
 31" Height, Adjustable to 45"
 31½" Tabletop Length
 12" Tabletop Width
 16" Overall Width
 3" Casters




INJECTIONSTAND

MODEL # IVM1925-01
 31½" Height, Adjustable to 50½"
 11½" Pad Length
 7" Pad Width
 21" Base Width
 2" Casters



MRI SOLID OAK COIL CARTS

Non-Magnetic Solid Oak Constructed
 5" Casters

NO ASSEMBLY REQUIRED 
 Cart Does Not Come With Cover or Handles
 Burgundy Cover or Handle Optional



**OPTIONAL CART COVER
 NOT INCLUDED
 ADJUSTABLE
 SHELVES**



PICTURED WITH OPTIONAL
 HANDLES & COVER

OPTIONAL HANDLES MODEL # FRM2028-01



Non-ferromagnetic Solid Oak Constructed
NO ASSEMBLY REQUIRED
 Cart Does Not Come With Cover or Handles
 Burgundy Cover or Handle Optional



**NON-FERROMAGNETIC
 CART CASTERS**

**MR CONDITIONAL
 TO 7 TESLA**

MRI NON-MAGNETIC	MRI NON-FERROMAGNETIC	SIZE	# OF SHELVES	SHIP WEIGHT	SHIP D X H X L
FRM2000-01	FRM2500-01	2'D x 4'H x 2'L	3 Shelves	145 lbs	4'D x 5'H x 4'L
FRM2001-01	FRM2501-01	2'D x 4'H x 3'L	3 Shelves	160 lbs	4'D x 5'H x 4'L
FRM2002-01	FRM2502-01	2'D x 4'H x 4'L	3 Shelves	185 lbs	4'D x 5'H x 4'L
FRM2003-01	FRM2503-01	2'D x 6'H x 2'L	4 Shelves	185 lbs	4'D x 7'H x 4'L
FRM2004-01	FRM2504-01	2'D x 6'H x 3'L	4 Shelves	185 lbs	4'D x 7'H x 4'L
FRM2005-01	FRM2505-01	2'D x 6'H x 4'L	4 Shelves	200 lbs	4'D x 7'H x 4'L
FRM2011-01	FRM2511-01	2'D x 3'H x 2'L	2 Shelves	130 lbs	4'D x 4'H x 4'L
FRM2012-01	FRM2512-01	2'D x 3'H x 3'L	2 Shelves	145 lbs	4'D x 4'H x 4'L
FRM2013-01	FRM2513-01	2'D x 3'H x 4'L	2 Shelves	160 lbs	4'D x 4'H x 4'L

MRI SOLID OAK SHELVING



MRI solid oak shelving comes in a variety of sizes.
All sizes have 6 legs, with or without casters.
The shelves are 2' Deep, 3' to 6' in Height, and 6' or 8' in Length.
The Top is considered as one of the shelves.
No Assembly Required
5" Casters

**SHELVING IS CUSTOM MADE AND NON-RETURNABLE
CALL FOR SHIPPING INFORMATION AND PRICING**

MODEL #	SIZE	# OF SHELVES	WITH CASTERS	SHIP WEIGHT	SHIP D X H X L
FRM2020-01	2'D x 3'H x 6'L	3 Shelves	Yes	145 lbs	4'D x 5'H x 4'L
FRM2021-01	2'D x 4'H x 6'L	3 Shelves	Yes	160 lbs	4'D x 5'H x 4'L
FRM2022-01	2'D x 6'H x 6'L	3 Shelves	Yes	185 lbs	4'D x 5'H x 4'L
FRM2023-01	2'D x 3'H x 8'L	3 Shelves	Yes	185 lbs	4'D x 7'H x 4'L
FRM2024-01	2'D x 4'H x 8'L	4 Shelves	Yes	185 lbs	4'D x 7'H x 4'L
FRM2025-01	2'D x 6'H x 8'L	4 Shelves	Yes	200 lbs	4'D x 7'H x 4'L
FRM2120-01	2'D x 3'H x 6'L	3 Shelves	No	200 lbs	4'D x 9'H x 4'L
FRM2121-01	2'D x 4'H x 6'L	3 Shelves	No	215 lbs	4'D x 9'H x 4'L
FRM2122-01	2'D x 6'H x 6'L	3 Shelves	No	230 lbs	4'D x 9'H x 4'L
FRM2123-01	2'D x 3'H x 8'L	3 Shelves	No	130 lbs	4'D x 4'H x 4'L
FRM2124-01	2'D x 4'H x 8'L	4 Shelves	No	145 lbs	4'D x 4'H x 4'L
FRM2125-01	2'D x 6'H x 8'L	4 Shelves	No	160 lbs	4'D x 4'H x 4'L



FRM-2024-01



MRITABLES 



The **MRi furniture** line is built completely with non-magnetic supplies. The furniture is made to go directly into the MR environment. Including: Chair, Sofa, End Table, and Coffee Table (not pictured). Choose from 4 fabric choices. Items are built to order, please inquire on lead times. Custom furniture available upon request.

MODEL #	DESCRIPTION	DIMENSIONS	SHIP WEIGHT	SHIP L X W X H
FRM1002-01	End Table	24" W x 24" D 24" High	305 lbs	25" x 25" x 25"
FRM1003-01	Coffee Table	24" W x 48" L 24" High	305 lbs	30" x 58" x 34"



MRIPEDIATRIC CHAIR & TABLE 

**PURCHASE 1 PEDIATRIC TABLE AND 2 PEDIATRIC CHAIRS
 SAVE \$50.00**



MODEL #	DESCRIPTION	DIMENSIONS	SHIP WEIGHT	SHIP L X W X H
FRM1005-01	Pediatric Chair	16" Wide x 14" Deep Seat Height: 10" No Arms	48 lbs	20" x 20" x 20"
FRM1006-01	Pediatric Table	21" Long x 25 1/2" Wide x 18" High	94 lbs	25" x 25" x 25"


MRI FURNITURE **CHAIRS & SOFA** 



FABRIC CHOICES




NAVY BLUE



BURGUNDY



FOREST GREEN



GREY



MODEL #	DESCRIPTION	DIMENSIONS	SHIP WEIGHT	SHIP L X W X H
FRM1000-01	Sofa	53" W x 24 1/2" D Seat Height: 20" Arm Height: 28"	305 lbs	30" x 58" x 34"
FRM1001-01	Chair	29" W x 24 1/4" D Seat Height: 20" Arm Height: 28"	305 lbs	36" x 36" x 36"
FRM1004-01	Bariatric Chair	31" W x 24 1/4" D Seat Height: 20" Arm Height: 28" 500 lb Weight Cap.	305 lbs	36" x 36" x 36"
FRM1005-01	Pediatric Chair	16" Wide x 14" Deep Seat Height: 10" No Arms	48 lbs	20" x 20" x 20"

MRI TABLES/CARTS 



INM1507-01



INM1519-01



INM1501-01



INM1513-01

WARRANTY

5 Years on Frame, 1 Year on Moving Parts

Constructed of Stainless Steel

2" Dual Wheel Casters

Assembly Required

MODEL #	SIZE	W/RAILS	# OF SHELVES	SHIP WEIGHT	SHIP L X W X H
INM1507-01	18" x 24"	No	1 Top	13 lbs	20" x 5" x 26"
INM1519-01	18" x 24"	No	2	23 lbs	20" x 5" x 26"
INM1501-01	18" x 24"	Yes	1 Top	18 lbs	20" x 5" x 26"
INM1513-01	18" x 24"	Yes	2	28 lbs	20" x 5" x 26"



ALL TABLES SHIP WITH 2ND BOX CONTAINING 4 LEGS AND HARDWARE.

MRI POLYMER UTILITY CART WITH 3 SHELVES AND 3" NON MAGNETIC CASTERS 



Enclosed end panels
Convenient handle on each end
Conditional to 3T

Specifications:

Length 34" (86.36 cm)

Width 18" (45.72 cm)

Height 37" (93.98 cm)

WEIGHT CAPACITY IS 150 LBS (68.03 KG)

MODEL # MTM1990-01

Ship Weight: 80 lbs

Ship Dimensions: 40" x 48" x 60"

6 DRAWER STAINLESS STEEL CART 



Dimensions:

27" Wide
42" Tall
16½" Deep
5" Casters



Drawer Dimensions:

11" x 21½" Drawers
Two 3" Deep Drawers
Three 4" Deep Drawers
One 8" Deep Drawer
One 11" x 13½" Table Drawer



MODEL # TAM1460-01

Ship Weight: 150 lbs

Ship Dimensions: 45" x 40" x 24"



PVC SPECIALTY CART W/4 SHELVES


MODEL # PVM1458-01

4 Shelves (Size 20x25)
8 Pull Out Tubs (Size 20"x15"x5")
Mesh or Solid Vinyl Cover with color choice.

All carts come with solid vinyl fabric, mesh can be substituted.
Carts shipped with Royal Blue vinyl unless color specified.

Height: 12" Between Shelves
Length: 48" L without handles X 65" H X 21.5" W
50.5" L with handles X 65" H X 21.5" W

Available Colors



PVC LINEN/MULTI-USE CARTS

FEATURES:

MR COMPATIBLE

Totally enclosed carts help prevent dust build-up on expensive coils.

All carts come with solid vinyl fabric, mesh can be substituted.
Carts are shipped with Royal Blue vinyl, unless color specified.

Velcro closures

Caster sockets are double wall reinforced with schedule 40 one inch PVC
Can upgrade to 5" casters for additional charge

ALL MULTI-USE CARTS ARE SHOWN WITH THE FRONT COVER OPEN.

MODEL	# OF SHELVES	SHELF SIZE	HEIGHT BETWEEN SHELVES	HEIGHT WITH CASTERS	WEIGHT CAPACITY PER SHELF
PVM1430-01	3	20" x 25"	12"	33"	75 lb
PVM1431-01	2	25" x 20"	12"	37"	75 lb
PVM1433-01	3	25" x 20"	12"	51 1/4"	75 lb
PVM1434-01	4	25" x 20"	12"	65 1/2"	75 lb
PVM1436-01	2	32" x 20"	12"	37"	100 lb
PVM1438-01	3	32" x 20"	12"	51 1/4"	100 lb
PVM1439-01	4	32" x 20"	12"	65 1/2"	100 lb
PVM1442-01	2	45" x 20"	16"	46"	125 lb
PVM1443-01	3	45" x 20"	16"	65"	125 lb
PVM1444-01	4	45" x 20"	16"	77"	125 lb
PVM1445-01	2	50" x 20"	16"	45 1/4"	150 lb
PVM1446-01	3	50" x 20"	16"	65"	150 lb
PVM1447-01	4	50" x 20"	16"	77"	150 lb





NON-FERROMAGNETIC LAUNDRY HAMPER



NEW FROM CREATIVE LIVING MEDICAL

Conceal and easily transport soiled linens and help support Infection Prevention. Solid Aluminum construction for bag rest, casters, and an optional foot pedal. It has been tested and certified to pass through the pre-screen and pre-entry door systems.

Hamper Bag Sold Separately.



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
HAA2300-01	Hamper w/ Lid	21 lbs	12" x 24" x 38"
HAA2301-01	Hamper w/ Lid and Foot Pedal	21 lbs	12" x 24" x 38"
HAA2503-01	Hamper Bag	5 lbs	11" x 11" x 5"

WARRANTY: 5 Years, 1 Year on Moving Parts

PVCHAMPERS








WARRANTY:

1 year on frame
6 months on casters and soft goods

All mesh hampers sent Royal Blue unless color specified. SEE COLORS on PG. 23
All leak proof hampers sent white unless color specified. No zipper on Leak Proof Bags



PVM1401-01


MODEL	# OF BINS	BAG TYPE	LENGTH	WIDTH	HEIGHT
PVM1400-01 	1	33 Gallon Mesh	18"	18"	37"
PVM1401-01 	1	14.46 Gallon Leak Proof	18"	18"	37"
PVM1407-01 	1	42 Gallon Mesh	26½"	22½"	35½"
PVM1402-01 	1	33 Gallon Mesh	23¾"	18¾"	38"
PVM1403-01 	2	22 Gallon Leak Proof	38¾"	18"	37"
PVM1404-01 	3	43 Gallon Mesh	58¾"	18"	37"
PVM1405-01 	3	67 Gallon Leak Proof	58¾"	22½"	35½"


STAINLESS STEEL HAMPERS


HAMPER BAG SOLD SEPARATELY


Height: 36"
Inside Diameter: 16 1/2"
Outside Diameter: 18"
Casters: 2" Dual Wheel

WARRANTY:
5 years on frame
1 year moving parts
90 Days on Soft Goods

MODEL # HAM2200-01
Hamper Only - No Lid 
Ship Weight: 75lbs
Ship Dimensions: 20" x 20" x 36"

MODEL # HAM2202-01 (PICTURED)
Hamper with Foot Operated Lid 
Ship Weight: 75lbs
Ship Dimensions: 20" x 20" x 36"

MODEL # HAM2201-01
Hamper with Lid 
Ship Weight: 75lbs
Ship Dimensions: 20" x 20" x 36"

MODEL # HAM2203-01
1 Hamper Bag 
Ship Weight: 3lbs
Ship Dimensions: 11" x 11" x 5"



FERROMAGNETIC DETECTION

Ferromagnetic detection systems have been demonstrated to be highly effective as a quality assurance tool, verifying the successful screening and identifying ferromagnetic objects which were not discovered by conventional screening methods. It is recommended that new facility construction anticipate the use of ferromagnetic detection screening and provide for installation of the devices in a location which facilitates use and throughput. The recommended use of ferromagnetic detection to verify the screening of patients before they pass through the controlled point of access into Zone III.

For utmost protection, deploy the screener as part of your comprehensive MRI ferromagnetic safety system; including our detection system guarding your magnet room doorway and ferrous-free patient transfer equipment.



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
MGM2510-01	FerrAlert II PLUS	45 lbs	8" x 13" x 85"

Ferromagnetic Detection System - MRI Entry System

Adverse MRI incidents are on the rise. While the number of incidents is small compared to the number of procedures performed, even one is too many for any clinic, hospital or research facility.

Our Ferromagnetic Detector protects the magnet room against moderate to large size ferromagnetic projectiles. This product significantly helps to protect the investment in the magnet and to avoid a catastrophic accident.



Operation

The **GREEN** LED array illuminated around the MRI entrance indicates that the detector is scanning and there are no ferromagnetic threats detected.

A **YELLOW** beacon light indicates that a threat is approaching the portal.

The **RED** LED array flashes the entire portal, indicating a major threat such as a ventilator or oxygen cylinder has approached the portal.

The greater the threat, the further away the warning is activated.

Ferrous objects not approaching the entrance to Zone IV do not trigger an audible warning.

Features

The highest sensitivity system available, detects even the smallest ferromagnetic objects.

Enables you to screen your patients without impacting workflow and existing team practices.

Becomes a simple step in your MRI safety routine.

Demonstrates to your patients just how serious you are about their safety and MRI image quality.

MRI NON-FERROMAGNETIC STRETCHER ⚠️
WITH FLIP DOWN SIDE RAILS
MODEL # STA5000-01

uspatentnumber
9,144,409 B1

FEATURES:

- NEW!** Aluminum rigid construction
- Drop down side rails
- 5" Ceramic Total Lock Casters
- Full surround rubber bumper
- Adjustable headrest
- 3" Pad

MR CONDITIONAL TO 7 TESLA



SPECIFICATIONS:

Height	29 1/2"
Height with Rails up	40"
Height of Rails	15"
Length	77 1/4"
Rail Length	39"
Width:	27"
Weight:	130 lbs.
Back Adj. Range	0-75 Degrees
Storage Tray:	1 1/2"D x 26"L x 17 1/2"W
Number of Positions:	11
WEIGHT CAPACITY:	500 LBS



MRI NON-FERROMAGNETIC STRETCHER ⚠️
WITH WRAP AROUND RAILS
MODEL # STA5300-01

WARRANTY:

- 5 Years on Frame
- 1 Year on Moving Parts
- 90 Days on Soft Goods

SPECIFICATIONS:

Height	30 1/4"
Height with Rails up	42"
Height of Rails	11 3/4"
Length	76"
Rail Length	61 7/8"
Width:	28 1/8"
Back Adj. Range	0-71.2 Degrees
Storage Tray:	1 5/8"D x 26"L x 17 1/2"W
Number of Positions:	10
WEIGHT CAPACITY:	500 LBS



NON-FERROMAGNETIC STRETCHER ACCESSORIES ⚠️ **ACCESSORIES**
MR CONDITIONAL TO 7 TESLA

I.V. POLE

MODEL # STA5100-01
Ship Weight: 12 lbs
Ship Dimensions: 8" x 8" x 36"

OXYGEN TANK HOLDER

MODEL # STA5101-01
Ship Weight: 12 lbs
Ship Dimensions: 8" x 8" x 36"

RAIL PADS (SOLD AS PAIR)

MODEL # STA1159-CL
Ship Weight: 25 lbs
Ship Dimensions: 14" x 14" x 24"



Specify Color



REPLACEMENT PAD

MODEL # STM1555-01-22 | BURGANDY
MODEL # STM1555-01-46 | BLACK
MODEL # STM1555-01-50 | GREEN
MODEL # STM1555-01-02 | DARK BLUE
Thickness: 3"
Width: 23"
Length: 73 1/2"
Ship Weight: 30 lbs
Ship Dimensions: 6" x 24" x 40"

REPLACEMENT CASTERS

MODEL # STM2353-01
Ship Weight: 21 lbs
Ship Dimensions: 14" x 6" x 4"
Sold Individually

NON-FERROMAGNETIC WHEELCHAIRS

FEATURES:

- Burgundy Upholstery
- Aluminum construction
- Flip-Up Arms
- Plastic Side Panel
- Swingaway Footrest
- WEIGHT CAPACITY: 300 LBS**

WARRANTY:

- 5 Years on Frame & Crossbraces
- 1 Year on Moving Parts
- 90 Days on Soft Goods

Ship Weight: 76 lbs

Ship Dimensions: 12" x 34" x 38"

SPECIFICATIONS:

- Seat Width: 20"
- Overall Width: 28"
- Overall Depth: 36" w/footrest
- Seat Depth: 16"
- Overall Height: 36"
- Back Height: 19"
- Front Casters: 8" x 1" solid
- Rear Wheel: 24" solid
- Wheel Locks: push to lock
- Weight: 31 lbs.
- with footrests
- 300 LB WEIGHT CAPACITY**

SPECIFICATIONS:

- Seat Width: 22"
- Overall Width: 30"
- Overall Depth: 43" w/footrest
- Seat Depth: 18"
- Overall Height: 37"
- Back Height: 19"
- Front Casters: 8" x 1" solid
- Rear Wheel: 24" solid
- Wheel Locks: push to lock
- Weight: 33 lbs.
- with footrests
- 300 LB WEIGHT CAPACITY**

SPECIFICATIONS:

- Seat Width: 24"
- Overall Width: 31½"
- Overall Depth: 44" w/footrest
- Seat Depth: 20"
- Overall Height: 36"
- Back Height: 19"
- Front Casters: 8" x 1" solid
- Rear Wheel: 24" solid
- Wheel Locks: push to lock
- Weight: 43½ lbs.
- with footrests
- 500 LB WEIGHT CAPACITY**

20" WIDE ALUMINUM WHEELCHAIR

- Detachable Footrests
- BURGUNDY MODEL # WCA2025SD27**
- GREEN MODEL # WCA2025SD50**
- DARK BLUE MODEL # WCA2025SD02**
- Padded Arms
- MR CONDITIONAL TO 7 TESLA**



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



Back of upholstery

22" WIDE ALUMINUM WHEELCHAIR

- Detachable Footrests
- BURGUNDY MODEL # WCA2226SD27**
- GREEN MODEL # WCA2226SD50**
- DARK BLUE MODEL # WCA2226SD02**
- Padded Arms
- MR CONDITIONAL TO 7 TESLA**



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

24" WIDE ALUMINUM WHEELCHAIR

- Detachable Footrests
- BURGUNDY MODEL # WCA2425SD27**
- GREEN MODEL # WCA2425SD50**
- DARK BLUE MODEL # WCA2425SD02**
- Padded Arms
- MR CONDITIONAL TO 7 TESLA**

NON-FERROMAGNETIC WHEELCHAIR ACCESSORIES



WARRANTY:
5 Year Warranty



ACM1503-01

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
ACM1600-01	Aluminum I.V. Pole	12 lbs	8" x 8" x 36"
ACM1601-01	Aluminum I.V. & Oxygen Tank Holder	12 lbs	8" x 8" x 36"
ACM1602-01	Aluminum Oxygen Tank Holder	10 lbs	10" x 10" x 19"
ACM1503-01	Stencil Company Name or Other Information	5 lbs	4" x 12" x 19"

NON-FERROMAGNETIC LAUNDRY HAMPER



WARRANTY:
5 Years on Frame
1 Year on Moving Parts

BAG SOLD SEPARATELY

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
HAA2300-01	Hamper w/ Lid	21 lbs	12" x 24" x 38"
HAA2301-01	Hamper w/ Lid and Foot Pedal	21 lbs	12" x 24" x 38"
HAA2503-01	Hamper Bag	5 lbs	11" x 11" x 5"

NON-FERROMAGNETIC MRI STEP STOOLS



500 LB CAPACITY
5 YEAR WARRANTY EXCLUDING TIPS



MODEL # FSA3000-01

Step Stool without Handrail
16" x 12" Base

Ship Weight: 18 lbs


Ship Dimensions: 14" x 14" x 18"

MODEL # FSA3001-01

Step Stool with 41" Handrail
16" x 12" Base

Ship Weight: 48 lbs

Ship Dimensions: 15" x 15" x 42"



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.

NON-FERROMAGNETIC MRI FOLDING WALKERS



MRI JUNIOR FOLDING WALKER

MODEL # WAA1001-01

Ship Weight: 10 lbs

Ship Dimensions: 24" x 5" x 24"

MRI ADULT FOLDING WALKER

MODEL # WAA1000-01

Ship Weight: 10 lbs

Ship Dimensions: 36" x 5" x 24"

MRI HEAVY DUTY / EXTRA WIDE FOLDING WALKER

MODEL # WAA1002-01

Ship Weight: 10 lbs

Ship Dimensions: 36" x 5" x 24"




WARRANTY:
5 Year Warranty
Excluding Tips on Walkers



Both the adult and heavy duty walker comply with the essential requirements of the relevant European health, safety and environmental protection legislation.



NON-FERROMAGNETIC IV POLES / O² HOLDERS



MRI NON-FERROMAGNETIC

IV POLES WITH BASE

Adjustable Base
Adjustable Height: 56" to 105"

MODEL # IVA3000-01

2 IV Hooks, Handle, 5
Non-Locking Casters

MODEL # IVM1311-01

2 IV Hooks, Handle, 5
Locking Casters

MODEL # IVA3001-01

2 IV Hooks, Handle, 5
Non-Locking Casters

MODEL # IVM1312-01

4 IV Hooks, Handle,
5 Locking Casters

MODEL # IVM1350-01

Steering Handle

MRI NON-FERROMAGNETIC

SINGLE O² TANK HOLDERS WITH BASE

Adjustable Base w/ 5 Casters
Adjustable Height: 56" to 105"

MODEL # OXA3000-01

No IV Hooks, Accessory Hook, Handle

MODEL # OXA3001-01

2 IV Hooks, Accessory Hook, Handle

MODEL # OXA3002-01

4 IV Hooks, Accessory Hook, Handle

MODEL # IVM1350-01

Steering Handle

WARRANTY:

5 years on frame
1 Year on Moving Parts

MRI NON-FERROMAGNETIC

DUAL O² TANK HOLDERS WITH BASE

Adjustable Base w/ 5 Casters
Adjustable Height: 56" to 105"

MODEL # OXA3003-01

No IV Hooks, Accessory Hook, Handle

MODEL # OXA3004-01

2 IV Hooks, Accessory Hook, Handle

MODEL # OXA3005-01

4 IV Hooks, Accessory Hook, Handle

MODEL # IVM1350-01

Steering Handle

CALL FOR SHIPPING DETAILS

DIMENSIONS:
Size: 9" x 13"



SET OF 4 PLASTIC **CLIPBOARDS**

MODEL #CLM1012-01
1 of each color: Green, Violet, Royal Blue, Orange
Clipboard Sheet Size: 8 1/2" W x 11" L
Clip Capacity: 1/2" of 20Lb paper with built in pen holder.

SET OF 4
Ship Weight: 6 lbs
Ship Dimensions: 19" x 12" x 4"

BREAST CANCER AWARENESS **CLIPBOARDS**

DIMENSIONS:
Size: 9 5/8" x 13 3/8"



MODEL # CLM1763-01
Clipboard Sheet Size: 8 1/2" Wide x 11" Long
Clip Capacity: 1/2" of 20 lb paper with built in pen holder
Pink Clipboard
Pink Ribbon Breast Cancer Awareness along bottom
Product was tested by a facility in a 3T magnet and was said to have very little torque right in the bore of the magnet.

SET OF 4
Ship Weight: 6 lbs
Ship Dimensions: 19" x 12" x 4"



NON-MAGNETIC **PAPER CLIPS**

Comes in assorted colors

MODEL #	SIZE	QUANTITY
MTM1702-01	Small	1000
MTM1703-01	Medium	500
MTM1701-01	Large	200
MTM1700-01	Extra Large	50



MRI COMPLETED **LABELS**

MODEL # MTM1460-01
Label Size: 1.875" x .875"
Great for use on patient folders.
Lime Green label with Black Lettering.

500 PER ROLL
Ship Weight: 4 lbs
Ship Dimensions: 7" x 4" x 4"



PILLOW **POSITIONERS**

Relieves stress on the shoulders during diagnostic MRI and CT procedures.
25" Long x 16" Wide x 8 1/4" High

MODEL # FOM1009-01
Flat Bottom
Ship Weight: 42lbs
Ship Dimensions: 20" x 20" x 20"

MODEL # FOM1010-01
Contoured Bottom
Ship Weight: 42lbs
Ship Dimensions: 20" x 20" x 20"

UTILITY RACK

MODEL # CRM1300-01

Stainless Steel

Height: 66 1/2" Tall

Ship Weight: 23lbs

Ship Dimensions: 8" x 8" x 67"



NON-FERROMAGNETIC SAND BAG POSITIONERS

Custom upholstered sand bag positioners with supportive handles.

MR SAFE TO 7 TESLA



MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
SBG1003-01	3 lb. Sand Bag Positioner, 8" W x 8" L	3 Lbs	4 Lbs	11" x 11" x 5"
SBG1005-01	5 lb. Sand Bag Positioner, 5.5" W x 15" L	5 Lbs	6 Lbs	11" x 11" x 5"
SBG1007-01	7 lb. Sand Bag Positioner, 9" W x 12" W	7 Lbs	8 Lbs	19" x 12" x 4"
SBG1010-01	10 lb. Sand Bag Positioner, 9" W x 14" L	10 Lbs	11 Lbs	19" x 12" x 4"
SBG1012-01	12.5 lb. Sand Bag Positioner, 12" W x 12"	12.5 Lbs	13.5 Lbs	19" x 12" x 4"
SBG1013-01	SET OF 5, 1 ea. of 3, 5, 7, 10, and 12.5 LB. Sand Bag Positioners	37.5 Lbs	38.5 Lbs	14" x 14" x 18"

SAND BAG POSITIONERS

Black outer bag is naghayde reinforced naghayde, and dielectrically sealed to prevent leakage or dust. Inner bag is electrically sealed for added protection.

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
FOM1200-01	7" W x 10" L, 1 Each	3 Lbs	5 Lbs	11" x 11" x 5"
FOM1201-01	6" W x 18", 1 Each	7 Lbs	9 Lbs	11" x 11" x 5"
FOM1203-01	6" W x 17" L, 1 Each	12 Lbs	14 Lbs	19" x 12" x 4"
FOM1204-01	11" W x 11" L, 1 Each	12 Lbs	14 Lbs	19" x 12" x 4"
FOM1205-01	Set of 4 Sand bags	34 Lbs	37 Lbs	18" x 10" x 8"

SAND BAG POSITIONERS

Colored outer bag is nylon reinforced vinyl, and dielectrically sealed to prevent leakage or dust. Inner bag is electrically sealed for added protection. Covers are nylon reinforced, color coded vinyl. Complete with convenient handle.

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
FOM1300-01	Yellow Sand Bag, 1 Each	3 Lbs	4 Lbs	11" x 11" x 5"
FOM1301-01	Orange Sand Bag, 1 Each	5 Lbs	6 Lbs	11" x 11" x 5"
FOM1302-01	Red Sand Bag, 1 Each	7 Lbs	8 Lbs	19" x 12" x 4"
FOM1303-01	Blue Sand Bag, 1 Each	10 Lbs	11 Lbs	19" x 12" x 4"
FOM1326-01	Set of 4 Sand bags	29 Lbs	30 Lbs	18" x 10" x 8"



This product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.



SWIRLEARPLUGS MR

Fun colored to spark up the hearing protection program.
Different colors in every package.
Tapered for easy insertion in any size ear canal.
Non-irritating, smooth surface for added comfort.

NOISE REDUCTION RATING 32

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
RAM1009-01	Without Cord, 1000 Count/500 Pair	3 Lbs	8 Lbs	18" x10" x 8"
RAM4006-01	With Cord, 200 Count/100 Pair	3 Lbs	8 Lbs	18" x10" x 8"



CONTAINERS ARE NON-REFILLABLE

PLUGDISPENSER MR

Clicking mechanism for easier use and dispensing.
Chute delivers plugs right into the palm of your hand.
Affordable enough to be placed near every work area.
Eliminates the mess of individual paper/polybag packaging.
Each dispenser is sealed and completely recyclable.
Available in 250 or 500 count per plug dispenser with either Swirl Plugs or Pure Fit foam earplugs.
Mounting bracket with hardware included with every plug dispenser.

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
RAM4009-01	Swirl, 250 Count/125 Pair	2 Lbs	8 Lbs	18" x10" x 8"
RAM4010-01	Swirl, 500 Count/250 Pair	3 Lbs	8 Lbs	18" x10" x 8"
RAM4011-01	Pure Fit, 250 Count/125 Pair	2 Lbs	8 Lbs	18" x10" x 8"
RAM4012-01	Pure Fit, 500 Count/250 Pair	3 Lbs	8 Lbs	18" x10" x 8"

DARK BLUE EYE SHADE **SLEEP MASK** MR

MODEL # **RAM6000-01 NYLON 100 PER CASE**
MODEL # **RAM6001-01 NYLON 1000 PER CASE**

MRI safe eyeshade sleep mask. Made from satin soft polyester, this is the same mask provided by major airlines and hotels. Available in dark blue only (inside is black).
2 elastic straps.

Ship Weight: 5 lbs

Ship Dimensions: 8" x 8" x 12"



NON MAGNETIC **PRISM GLASSES** MR

MODEL # **RAM2013-01**

Features:

Fits over Prescription Glasses.

Do not reverse the image so it is possible to see a mural or projection screen TV.

Prism angle remains fixed in either position so no viewing angle adjustment is required.

Free of any metal components and may be worn in all types of Diagnostic Imaging Scans.

Similar to wearing a conventional pair of glasses so the patient will easily adapt to wearing them.

Ship Weight: 2 lbs

Ship Dimensions: 4" x 4" x 7"



NEONATAL **NOISE GUARDS** MR

A must have for premie and newborn patients, Neonatal Noise Guards protect these little one's sensitive ears in the MRI environment, and provide peace and quiet in the NICU. Neonatal Noise Guards reduce the surrounding environment noise level by at least 7db, with a reduction in sound pressure of more than 50%.

Disposable and latex free.

BOX INCLUDES 6 PAIRS

MODEL # **RAM1016-01**

Neonatal Noise Guards

6 Pair Per Box

Ship Weight: 3 lbs

Ship Dimensions: 6" x 6" x 6"



NEONATAL NOISE GUARD **DISPENSER** MR

MODEL # **RAM1030-01**

Makes it easy to store and quickly access Neonatal Minimuffs or any neonatal ear plugs in a similar sized package.

Specifications:

Gravity Fed Dispenser

Slotted Thumb Opening

Holds Up To 12 Packages (Pairs)

Wall Mount With Brass Screws (Included)

3.75"L X 6.75"W x 6"H

Ship Weight: 2 lbs

Ship Dimensions: 7" x 8" x 10"



SMALL OR LARGE **SANITARY HEADSET** MR **COVERS** WHITE OR BLACK

MRI Non-Magnetic Sanitary Headset Covers prolong the life of headphones by keeping them clean and protect earmuffs from sweaty and unpleasant smells. Our headset cover will absorb oils and is a good hygienic practice in general. The elastic and stretchable material is ideal for a perfect fit. Prevent flakey headphone ear foam pad and earpieces on your headphones from breaking down, thus preventing the spread of infection. Our headset covers have no sound interference.

1000 PER CASE

MODEL # **RAM1007-01** 2 - 3" Diameter White

MODEL # **RAM1006-01** 5 - 6" Diameter White

MODEL # **RAM1058-01** 5 - 6" Diameter Black

Ship Weight: 5 lbs

Ship Dimensions: 8" x 8" x 12"



MRI NON-MAGNETIC **EAR MUFFS** MR

MODEL # **RAM1002-01** MRI Slimline Headset (29db) **3 YR WARRANTY**

MODEL # **RAM1057-01** Ultra-33 Passive Muffs (33db) **5 YR WARRANTY**

MODEL # **RAM1013-01** MRI Slimline Ear Muffs (25db) **90 DAY WARRANTY**

MODEL # **RAM1027-01** Replacement Headband

Magnetic Resonance Imaging (MRI) machines are very noisy (and sometimes even scary) places to spend time. Fortunately, most MRI labs now offer patients the option to wear ear muffs while inside the machine, where the noise is loudest.



RAM1002-01



RAM1013-01



RAM1057-01




RAM1027-01

MRI STEREO HEADPHONE ADAPTER

MODEL #RAM1081-01

MRI Stereo Headphone Adapter is designed for use MRI scanners that have their own existing patient sound system and a headphone jack that looks like a small 90-degree angled plug. The single mono plastic plug style headphone port is common to most of the Siemens MRI scanners. The single tube port audio connection provides mono audio for the MRI headphones.



MRI-SAFE STETHOSCOPE STYLE PNEUMATIC HEADPHONE MR

MODEL #BPM1066-01

MRI-Safe Stethoscope Style Pneumatic Headphone is ideal for MRI sound systems. It is completely made from plastic and foam so it will not interfere with MRI scans. The Stethoscope style MRI headphone will fit into most MRI coils, including head coils. Sound travels through plastic tubing to the left and right headphone earpieces. The ear tips may be replaced after each patient so the headphones remain clean and fresh for repetitive use.



MRI-SAFE HEARING PROTECTIVE EARMUFFS MR

MODEL #RAM1077-01

MRI-Safe Hearing Protective Earmuffs are very effective in reducing the high levels of noise produced by MRI systems. They help to protect MRI patients from hearing damage. MRI-Safe earmuffs have an adjustable headband to fit adults, teens and children. Each padded earpiece contains 2 layers of cellular foam to absorb additional noise. Earmuff is certified MRI SAFE totally free of any metal components and safe for all MRI scanners. Noise reduction rating of 20 dB.



IN-EAR MRI HEADPHONES MR

MRIaudio's newly designed iN-Ear Headphones are the ultimate solution for providing MRI patients with music, communication, and hearing protection. Their low profile design allows them to fit into every head coil and disposable foam ear tips prevent the spread of infection.

Adjustable

The optional/removable headband adjusts to fit any size head, while the swiveling ear tip inserts fit into every ear perfectly!

Versatility

The standardized pneumatic tube connector means that they will work with most MRI patient stereos available.

iN-Ear Headphones are **MR Safe**

Comes with a 1 year warranty.

Specifications:

- Noise reduction rating of 29db
- Clear, consistent sound quality
- Compressible, high-density foam eartips
- 0.5" headband profile fits in most MRI head coils
- Adjustable headband fits most heads



RAM3700-01



RAM3610-01
RAM3611-01

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
RAM3700-01	iN-Ear MRI Headphones	2 Lbs	4 Lbs	11" x 11" x 5"
RAM3611-01	Replacement Ear Tips/250 Pair	3 Lbs	8 Lbs	13" x 8" x 8"
RAM3610-01	Replacement Ear Tips/500 Pair	4 Lbs	8 Lbs	18" x 10" x 8"



GRE2000-01

INVISILUBE HEAVY DUTY CLEAR AEROSOL GREASE

Invisilube is designed for those specific applications where a clear lubricant is necessary. This clear grease spray works well with wheelchair bearings, journals and other parts that require lube.

MODEL #	DESC.	DESIGNATION
GRE2000-01	Invisilube Grease	

GentleGuard **Vinyl** 
Powder Free Gloves

Packed: 100 gloves per Easy Open Dispenser Box

MODEL#	DESCRIPTION
GLM2000-01-XS	Extra Small
GLM2000-01-SM	Small
GLM2000-01-MD	Medium
GLM2000-01-LG	Large
GLM2000-01-XL	Extra Large



GentleGuard **Latex** 
Powder Free Gloves

Packed: 100 gloves per Easy Open Dispenser Box

MODEL#	DESC.
GLM2001-01-XS	Extra Small
GLM2001-01-SM	Small
GLM2001-01-MD	Medium
GLM2001-01-LG	Large
GLM2001-01-XL	Extra Large



KN95 Mask **FDA & CE Certified** 

- MR Conditional at 3T or less
- With no artifact and it does not set off the solo screener

FREE SHIPPING ON ALL KN95 MASKS

MODEL#	DESC.
FS3000X-010	Case of 10
FS3000X-050	Case of 50
FS3000X-100	Case of 100
FS3000X-500	Case of 500



MRI WARNING SIGNS 

12" X 18" SIGNS IN SPANISH AND ENGLISH



MTM1350-01




MTM1351-01



MTM1358-01



MTM1359-01



SIGN DETAIL OF: MTM1350-01 & MTM1351-01



MTM1366-01



MTM1367-01

4" X 8" STICKERS AND SIGNS IN SPANISH AND ENGLISH

CUSTOM SIGNS AVAILABLE - CALL FOR MORE DETAILS


New codes are requiring facilities to use **Reflective** Warning Stickers and Signs in case of fire or power outage*.

MODEL #	SIZE	LANGUAGE	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
MTM1350-01	12" x 18"	English	Plastic Sign - Items Not Allowed	6 lbs	19" x 12" x 4"
MTM1442-01	12" x 18"	English	Reflective Plastic Sign - Items Not Allowed	6 lbs	19" x 12" x 4"
MTM1351-01	12" x 18"	Spanish	Plastic Sign - Items Not Allowed	6 lbs	19" x 12" x 4"
MTM1443-01	12" x 18"	Spanish	Reflective Plastic Sign - Items Not Allowed	6 lbs	19" x 12" x 4"
MTM1358-01	12" x 18"	English	Plastic Sign - No Unauthorized Entry	6 lbs	19" x 12" x 4"
MTM1444-01	12" x 18"	English	Reflective Plastic Sign - No Unauthorized Entry	6 lbs	19" x 12" x 4"
MTM1359-01	12" x 18"	Spanish	Plastic Sign - No Unauthorized Entry	6 lbs	19" x 12" x 4"
MTM1445-01	12" x 18"	Spanish	Reflective Plastic Sign - No Unauthorized Entry	6 lbs	19" x 12" x 4"
MTM1366-01	4" x 8"	English	Plastic Sign - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1438-01	4" x 8"	English	Reflective Plastic Sign - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1367-01	4" x 8"	Spanish	Plastic Sign - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1439-01	4" x 8"	Spanish	Reflective Plastic Sign - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1368-01	4" x 8"	English	Sticker - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1440-01	4" x 8"	English	Reflective Sticker - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1369-01	4" x 8"	Spanish	Sticker - Magnet is Always On	2 lbs	4" x 4" x 10"
MTM1441-01	4" x 8"	Spanish	Reflective Sticker - Magnet is Always On	2 lbs	4" x 4" x 10"

*CODES ARE EFFECTIVE IN VARYING STATES, CONTACT YOUR LOCAL FIRE DEPT. FOR MORE INFORMATION.

MRI ZONE SIGNS & GENERAL DIAGRAM

Ship Weight: 3 lbs
Ship Dimensions: 4" x 12" x 19"



MTM1161-01

MODEL #	DESCRIPTION	DES.
MTM1152-01	10" x 14", Zone 1 Sign, 1 Each	
MTM1174-01	5" x 7", Zone 1 Sign, 1 Each	
MTM1250-01	10" x 14", Spanish Zone 1 Sign, 1 Each	
MTM1153-10	10" x 14", Zone 2 Sign, 1 Each	
MTM1175-01	5" x 7", Zone 2 Sign, 1 Each	
MTM1251-01	10" x 14", Spanish Zone 2 Sign, 1 Each	
MTM1154-01	10" x 14", Zone 3 Sign, 1 Each	
MTM1176-01	5" x 7", Zone 3 Sign, 1 Each	
MTM1252-01	10" x 14", Spanish Zone 3 Sign, 1 Each	
MTM1155-01	10" x 14", Zone 4 Sign, 1 Each	
MTM1177-01	5" x 7", Zone 4 Sign, 1 Each	
MTM1253-01	10" x 14", Spanish Zone 4 Sign, 1 Each	
MTM1160-01	10" x 14", All 4 Zone Signs, Set of 4	
MTM1178-01	5" x 7", All 4 Zone Signs, Set of 4	
MTM1260-01	10" x 14", All 4 Zone Signs in Spanish, Set of 4	
MTM1161-01	12" x 18", Zone Sign General Diagram, 1 Each	



MTM1152-01
MTM1174-01



MTM1250-01



MTM1153-01
MTM1175-01



MTM1251-01



MTM1154-01
MTM1176-01



MTM1252-01



MTM1155-01
MTM1177-01



MTM1253-01




MTM1160-01
MTM1178-01



MTM1260-01

MRI WARNING SIGN-COMBINATION #1

MODEL # MTM1136-01
MRI Warning Sign-Combination of MTM1350-01 & MTM1351-01.
For **Items Not Allowed in MRI Room**. Half English/Half Spanish 18 x 24, Plastic



SIGN DETAIL OF:
MTM1350-01 & MTM1351-01



MTM1136-01

MRI WARNING SIGN-COMBINATION #2

MODEL # MTM1137-01
MRI Warning Sign-Combination of MTM1358-01 & MTM1359-01.
People Not Allowed Unattended, Half English/ Half Spanish, 18 x 24, Plastic

MODEL # MTM1136-01 & MTM1137-01
Shipping Weight: 10 lbs, Shipping Dimensions: 4" x 18" x 24"



MTM1137-01

ASTM SAFETY STANDARDS FOR MR ENVIRONMENTS

The highly anticipated Standard Practice for Marking Medical Devices for Safety in the Magnetic Resonance Environment by American Society for Testing and Materials (ASTM) has arrived.

As recommended by ASTM Intl; the intent of this practice is to provide needed information about the safety of items in and near MR scanners using a compact and easily recognized set of symbols and terms. The terms MR safe and MR compatible as first defined in 1997 in the FDA draft guidance document, "A Primer on Medical Device Interactions with Magnetic Resonance Imaging Systems," have been used to describe the safety of devices in and near MR systems. However, there has been a great deal of confusion surrounding this terminology. An incorrect assumption is often made, and users often incorrectly assume that items labeled MR safe or MR compatible are safe or compatible for any MR environment. MR environments vary in terms of magnetic field strength and RF conditions. This new system eliminates much of that confusion.

ICON DESIGNATION Products **NOT** marked with a designation have been tested under the past methods and are currently being re-evaluated under the new ASTM Standards. Products will be updated continually.



MR Safe

Any item that poses no known hazards in all MR environments.

MR Conditional



Any item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use. Field conditions that define the specified MR environment include field strength, spatial gradient, dB/dt (time rate of change of the magnetic field), radio frequency (RF) fields, and specific absorption rate (SAR). Additional conditions, including specific configurations of the item, may be required.



MR Unsafe

Any item that is known to pose hazards in all MR environments.

DANGER! NOT MRI SAFE WARNING STICKERS

DIE CUT - PERMANENT ADHESIVE - LAMINATED



DANGER!
NOT
MRI Safe



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
MTM1430-01	Red - 5 Large, 5 Medium, 5 Small	2 Lbs	6" x 6" x 6"
MTM1431-01	Red - 6 Large, 9 Medium, 15 Small	2 Lbs	6" x 6" x 6"

SIZE	10 PACK MODEL #	50 PACK MODEL #	100 PACK MODEL #	SHIP WEIGHT	SHIP L X W X H
1"x3" Red	MTM1424-01	MTM1425-01	MTM1426-01	5 Lbs	6" x 6" x 6"
1½" x 6" Red	MTM1402-01	MTM1403-01	MTM1404-01	5 Lbs	6" x 6" x 6"
4" x 6" Red	MTM1427-01	MTM1428-01	MTM1429-01	5 Lbs	6" x 6" x 6"

FERROMAGNETIC DETECTION IDENTIFICATION

HANGING EQUIPMENT **SIGNS** EQUIPMENT **BADGES**



Identify your MR facility equipment with these large plastic placards. This visual aid ensures technicians to remove magnetic equipment before entering MR room.

MODEL # MTM1198-01

Plastic
White / Red / Black
8"H x 10"W

Ship Weight: 2 lbs
Ship Dims: 15" x 13" x 1"



MTM1950-01



MTM1951-01

Clearly identify your MR facility equipment. Non-Tearable Vinyl or Rigid Plastic Badges. Slot is .125 x 1.5 6"H x 4"W Pack of 25

MODEL #	DESCRIPTION	SHIP WEIGHT
MTM1950-01	Vinyl Tag, This Item Passes	2 Lbs
MTM1951-01	Vinyl Tag, This Item Does Not Pass	2 Lbs
MTM1952-01	Plastic Tag, This Item Passes	2 Lbs
MTM1953-01	Plastic Tag, This Item Does Not Pass	2 Lbs

MRSAFE STICKERS

Die Cut - Permanent Adhesive - Laminated

Ship Weight: 1 lbs

Ship Dimensions: 4" x 4" x 10"

MODEL #	DESCRIPTION
MTM4100-10	1½" x 2", 10 Pack
MTM4100-50	1½" x 2", 50 Pack
MTM4100-100	1½" x 2", 100 Pack
MTM4101-10	3½" x 4", 10 Pack
MTM4101-50	3½" x 4", 50 Pack
MTM4101-100	3½" x 4", 100 Pack



MRCONDITIONAL STICKERS

Die Cut - Permanent Adhesive - Laminated

Ship Weight: 1 lbs

Ship Dimensions: 4" x 4" x 10"

MODEL #	DESCRIPTION
MTM4102-10	1½" x 2", 10 Pack
MTM4102-50	1½" x 2", 50 Pack
MTM4102-100	1½" x 2", 100 Pack
MTM4103-10	3½" x 4", 10 Pack
MTM4103-50	3½" x 4", 50 Pack
MTM4103-100	3½" x 4", 100 Pack



MRNOT SAFE STICKERS

Die Cut - Permanent Adhesive - Laminated

Ship Weight: 1 lbs

Ship Dimensions: 4" x 4" x 10"

MODEL #	DESCRIPTION
MTM4104-10	1½" x 2", 10 Pack
MTM4104-50	1½" x 2", 50 Pack
MTM4104-100	1½" x 2", 100 Pack
MTM4105-10	3½" x 4", 10 Pack
MTM4105-50	3½" x 4", 50 Pack
MTM4105-100	3½" x 4", 100 Pack



MRMULTI-PACKS

Die Cut - Permanent Adhesive - Laminated

Ship Weight: 1 lbs

Ship Dimensions: 4" x 4" x 10"

MODEL #	DESCRIPTION
MTM4053-01	3 of Each Style, 18 Total
MTM4054-01	6 of Each Style, 36 Total
MTM4056-01	10 of Each Style, 60 Total
MTM4055-01	20 of Each Style, 120 Total
MTM8009-01	1/2" x 3/4", 50 of Each Style
MTM8011-01	Set of 25 1/2" x 3/4" Each Style



MRI FLOORSTICKERS



Do Not Cross

Magnet Always On

MODEL # MTM1929-01

MODEL # MTM1386-01

17" Diameter Floor Stickers

Ship Weight: 6 lbs

Ship Dimensions: 19" x 12" x 4"

"DANGER RESTRICTED ACCESS"

GAUSS LINESTICKER

8" Diameter

White / Red / Black



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
MTM2500-01	5 Gauss Line	6 lbs	19" x 12" x 4"
MTM2501-01	10 Gauss Line	6 lbs	19" x 12" x 4"
MTM2502-01	30 Gauss Line	6 lbs	19" x 12" x 4"
MTM2503-01	100 Gauss Line	6 lbs	19" x 12" x 4"

MRIZONE IV 17" FLOORSTICKER

MODEL # MTM1923-01

Reads: DANGER! MRI Zone IV, Screened MRI Patients under constant direct supervision of trained MRI personnel only
17" Diameter

White / Red / Black



CAUTION RESTRICTED TO SCREENED MRI PATIENTS



MR Black / Yellow, 4" x 11"

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
MTM1066-01	Sticker	3 Lbs	4 Lbs	18" x 10" x 8"
MTM1065-01	Plastic Sign	3 Lbs	4 Lbs	18" x 10" x 8"

DANGER MUST BE ACCOMPANIED BY TRAINED MRI PERSONNEL



MR Black, Red and White, 4" x 11"

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
MTM1068-01	Sticker	3 Lbs	4 Lbs	18" x 10" x 8"
MTM1067-01	Plastic Sign	3 Lbs	4 Lbs	18" x 10" x 8"

NOTICE APPROACHING MRI RESTRICTED AREA



MR Blue and White, 4" x 11"

MODEL #	DESCRIPTION	WEIGHT	SHIP WEIGHT	SHIP L X W X H
MTM1069-01	Sticker	3 Lbs	4 Lbs	18" x 10" x 8"
MTM1070-01	Plastic Sign	3 Lbs	4 Lbs	18" x 10" x 8"

FLOORSTICKER/CARPET WARNING SIGN "MAGNET ALWAYS ON"



MR Size: 3' x 5'

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
MTM1128-01	Sticker, Powerful Magnet Always on!	15 lbs	36" x 6" x 6"
MTM1129-01	Carpet, Restricted Area, MRI Zone IV!	15 lbs	36" x 6" x 6"
MTM1370-01	Carpet, Powerful Magnet Always on!	15 lbs	36" x 6" x 6"

Pressure Sensitive Vinyl, Self adhesive
Laminate is Flextrack, Non-skid

LIGHTED SIGNS



MR Size of Sign: 11"W x 9"T x 2 1/4"D
Mounting Bracket: 5 3/4" x 4 1/2" x 1 1/8" Tall


MODEL #	STYLE	LIGHT	MESSAGE	SHIP WEIGHT	SHIP L X W X H
MTM1331-01	One-Sided	Red	Magnet In Use	6 lbs	19" x 12" x 4"
MTM1333-01	One-Sided	Red	Magnet Always On	6 lbs	19" x 12" x 4"
MTM1337-01	Two-Sided	Red	Magnet In Use	6 lbs	19" x 12" x 4"
MTM1339-01	Two-Sided	Red	Magnet Always On	6 lbs	19" x 12" x 4"

FOR USE OUTSIDE MR ROOM

NOT MR SAFE
THIS PRODUCT IS MAGNETIC

I.V., VENTILATOR & EQUIPMENT STANDS 

Satinless Steel I.V. Stands //



Warranty:
2 Years on IV Poles
6 Months on Soft Goods
**Stands Do Not include
Cylinders or Ventilator**

**Model
IVM1307-01**
4 Hook
I.V. Pole

**Model
IVM1312-01**
4 Hook
I.V. Pole

Model # IVM1306-01
2 Hook I.V. Pole

Model # IVM1311-01
2 Hook I.V. Pole

MODEL # DESCRIPTION

OXM1701-01	Magellan & Tau Omni Ventilator Stand w/ Dual D/E Cylinder Holder & Top Tray, Aluminum Base, 1 1/4" Diameter Pole
OXM1708-01	Ventilator Stand w/ Dual D/E Cylinder Holders, Aluminum Base, 60" High, 1 1/4" Diameter Pole
IVM1306-01	2 Hook I.V. Pole, Aluminum Base w/ Non-Locking Casters, Adjustable Height 59" - 88"
IVM1307-01	4 Hook I.V. Pole, Aluminum Base w/ Non-Locking Casters, Adjustable Height 59" - 88"
IVM1311-01	2 Hook I.V. Pole, Aluminum Base w/ Locking Casters, Adjustable Height 59" - 88"
IVM1312-01	4 Hook I.V. Pole, Aluminum Base w/ Locking Casters, Adjustable Height 59" - 88"
OXM1761-01	Equipment Stand Only, 48" High, 1 1/4" Diameter Pole, Rubber Tips
OXM1762-01	Equipment Stand Only, 60" High, 1 1/4" Diameter Pole, Rubber Tips
OXM1759-01	Equipment Stand Only, 48" High, 1 1/4" Diameter Pole, Non-Locking Casters
OXM1760-01	Equipment Stand Only, 60" High, 1 1/4" Diameter Pole, Non-Locking Casters
OXM1775-01	Equipment Stand Only, 48" High, 1 1/4" Diameter Pole, Locking Casters
OXM1776-01	Equipment Stand Only, 60" High, 1 1/4" Diameter Pole, Non-Locking Casters

BOX 1

SHIPS IN 2 BOXES

BOX 2

EQUIPMENT	SHIP WEIGHT	SHIP L X W X H
Ventilator Stands	20 lbs	15" x 15" x 48"
60" Equipment Stand	23 lbs	9" x 9" x 61"
48" Equipment Stand	20 lbs	9" x 9" x 50"

EQUIPMENT	SHIP WEIGHT	SHIP L X W X H
Solid Stainless Base	13 lbs	37" x 6" x 25"
Aluminum Base	6 lbs	25" x 6" x 25"

OXYGEN CYLINDER CARTS & STANDS




MODEL # OXM8505-01
Wall Mount Cylinder Holder



MODEL # OXM8506-01
Cylinder Floor Stand



MODEL # OXM1700-01
Oxygen Cart
Stainless Steel



MODEL # OXM8500-01
Oxygen Cart
Aluminum



MODEL # OXM8501-01
Oxygen Cart
Dual Cylinders



MODEL # OXM8502-01
Oxygen Cart
4 Cylinders



MODEL # OXM8503-01
Oxygen Cart
6 Cylinders



MODEL # OXM8504-01
Deluxe Oxygen Cart
M/M60 Cylinder



MODEL # OXM8507-01
Oxygen Cart
Aluminum

MODEL #	DESCRIPTION	SPECIFICATIONS	SHIP WEIGHT	SHIP L X W X H
OXM8505-01	Wall Mount Cylinder Holder	Aluminum, Fits E Cylinder Tanks	3 lbs	6" x 8" x 11"
OXM8506-01	Free Standing Cylinder Holder	Aluminum Oxygen Stand - Holds all Size Cylinders, Free Standing with Plastic Securing Strap 14" x 14" x 12.25"	25 lbs	19" x 19" x 13"
OXM8500-01	Oxygen Cart	Aluminum, Holds D or E size tanks 38" H x 12" W x 10" L	23 lbs	13" x 11" x 31"
OXM8501-01	Dual Tank Cart	Aluminum, Holds Dual D or E Cylinders. 38.5" High	44 lbs	15" x 14" x 40"
OXM8502-01	Oxygen Cart	Aluminum Oxygen Cart for 4 D or E Cylinders, on Wheels 38.5" High x 13.75" Wide x 13.25" Long	47 lbs	15" x 15" x 40"
OXM8503-01	Oxygen Cart	Aluminum Oxygen Cart for 6 D or E Cylinders, on Wheels 40.25" High x 13.75" Wide x 16.75" Long	65 lbs	19" x 16" x 41"
OXM8504-01	Deluxe Oxygen Cart	Aluminum with Wheels, Holds M/M60 size tanks 43.5" High x 16" Wide x 15" Long	64 lbs	18" x 16" x 43"
OXM8507-01	Oxygen Cart	Aluminum with Wheels, Holds M (8") Size Tanks 38" High x 14" Wide x 13" Long	44 lbs	15" x 14" x 40"
OXM1700-01	Oxygen Cart	Stainless Steel Oxygen Cart for D or E Cylinders, Adjustable Handle 38½" to 44", 6" Plastic Wheels	28 lbs	14" x 12" x 25"


OXYGEN CYLINDERS & WRENCHES

MODEL #	DESCRIPTION	SPECIFICATIONS	SHIP WEIGHT	SHIP L X W X H
OXM1728-01	E Cylinder	25" High x 4 ¹ / ₂ " Dia. Uses Yoke Regulators	12 lbs	8" x 8" x 36"
OXM1729-01	M Cylinder	36 ¹ / ₄ " H x 8" Dia. Uses Nut & Nipple Regs	12 lbs	8" x 8" x 36"
OXM1757-01	Cylinder Wrench	Large and Small Cylinder	2 lbs	7" x 4" x 4"
OXM1731-01	Cylinder Wrench	Small Cylinders	2 lbs	7" x 4" x 4"




BLACK
STAINLESS STEEL
OXM1757-01




BLACK PLASTIC
OXM1731-01



OXM1728-01

OXM1729-01

DU-O-VAC PLUS SUCTION SYSTEM

0-350 mm hg of suction, 2 each 50 PSI power take offs, 500 ml collection canister with built in filter. Powered by 50 PSI gas source

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
OXM1732-01	DISS Wall Connection	3 lbs	6" x 6" x 14"
OXM1733-01	Ohmeda Wall Connection	3 lbs	6" x 6" x 14"
OXM1739-01	Chemtron Wall Connection	3 lbs	6" x 6" x 14"
OXM1734-01	Puritan Bennett Wall Connection	3 lbs	6" x 6" x 14"
OXM1738-01	Packs of 10. Polystyrene & K resin construction. Frosted writing strip with 100 ml increments up to 400 ml. Built in hydrophobic filter cleans 99% of contaminated gas also acts as a safety filter.	11 lbs	11" x 20" x 9"



REPLACEMENT
CANISTERS
OXM1738-01

DU-O-VAC PLUS CODE CART SUCTION SYSTEM

INCLUDES:

Du-O Vac Plus, Regulator, Cylinder Cart, E Cylinder Option, 0-15 lpm Flowmeter

Ship Weight: 36 lbs

Ship Dimensions: 13" x 13" x 41"

MODEL #	DESCRIPTION
OXM1735-01	Du-O-Vac Plus Code Cart w/ E Cylinder
OXM1736-01	Du-O-Vac Plus Code Cart w/o E Cylinder



DU-O-VAC PLUS WITH REGULATOR AND FLOWMETER SUCTION SYSTEM

MODEL # OXM1737-01

INCLUDES:

Main body, collection canister, O2 flowmeter, oxygen regulator.

(Does not include cylinder and stand)

Ship Weight: 15 lbs

Ship Dimensions: 9" x 11" x 21"





WIDE UNIVERSAL **ARMBOARD**

Accommodates both arm and supplies. Slippery bottom allows for easy insertion under patient's shoulder with minimal patient movement. Thin enough to slide without discomfort. Patient's weight holds it in place.

MODEL # FOM1404-01



ACCUFIT MRI **BREAST COIL DISPOSABLE DRAPE** STYLE TWO

100 PER CASE

MODEL # BSM1007-01



When used as a system, the draping components offer comfort and protection in a very user-friendly design. The coil draping can be utilized with or without the breast coil pad. This gives the patient more comfort and the technologist better accessibility to the breast.

COMPATIBLE WITH:

- Siemens 1.5T Avanto/Espree 7 Channel
- Siemens 1.5T Symphony 7 Channel
- Phillips 1.5T Achieva 7 Channel
- Phillips 3.0T SENSE
- May fit coils other than those listed

Ship Weight: 14 lbs

Ship Dimensions: 14" x 14" x 24"

- GE 8 Channel coil
- GE 1.5T Excite 7 Channel
- GE 3.0T 7 Channel Breast
- CP Siemens Breast Array
- Insight

ACCUFIT MRI SENTINELLE BREAST COIL DISPOSABLE DRAPE MR

These drapes will minimize the daily cleaning routine of coils as well as ease the mind of clientele. The drapes provide comfort to patients, while providing protection from contagious diseases such as MRSA. Provides complete coverage and protection of all Sentinelle Breast Coil Products.

TOTAL COVERAGE:

Added an extra 60" in length now covering the entire table which eliminates costly usage of sheets and provides quick procedure turn around.

25 PER CASE

MODEL # BSM1011-01

COMPATIBLE WITH:

Siemens MAGNETOM® Avanto

MAGNETOM Espree

MAGNETOM Symphony

GE Signa® HDx or HDi MRI

Toshiba 1.5T Vantage Atlas™

Vantage Titan™

May fit coils other than those listed

Ship Weight: 7 lbs

Ship Dimensions: 12" x 10" x 12"



ACCUFIT MRI HEADREST DISPOSABLE COVER MR

Keep your headset clean from make-up, perspiration, open wounds etc. A quick and easy way to maintain your equipment.

100 PER CASE

MODEL # BSM1003-01

Ship Weight: 8 lbs

Ship Dimensions: 14" x 10" x 9"



ACCUFIT MRI SENTINELLE MR LARGE HEADREST COVER

These drapes will minimize the daily cleaning routine of coils as well as ease the mind of clientele. The drapes provide comfort to patients, while providing protection from contagious diseases such as MRSA. Finally, your staff will be delighted with your decision to simplify clean up between patients, thus achieving a greater throughput.

100 PER CASE


MODEL # BSM1012-01

This product fits the large 5 inch high headrest utilized with the Sentinelle Breast Coil Table and stand alone coils.

Ship Weight: 8 lbs

Ship Dimensions: 12" x 8" x 11"





MRI BIOPSY TRAYS

SAFETY: The only MRI Non-Ferromagnetic Certified Biopsy Tray. Designed and certified within the guidelines of the new ASTM standards in effect so your facility no longer has to assume the liability for any safety accidents due to the items made with ferrous materials.

USAGE FLEXIBILITY WITH LESS WASTE: Based on customer feedback for increased ordering and procedure flexibility you'll note that we packaged the main biopsy components separately as well as in a MRI biopsy tray. Trays can also be used as a MRI localization tray.

UNIQUE DESIGN: The only MRI NON-Ferromagnetic certified disposable scalpel and MRI NON-Ferromagnetic certified needle.

COST EFFECTIVE: This product is a first of its kind. Developed by technologists who perform the procedures. It is cost effective, practical and safe to use. The trays are completely sterile and disposable.

MR CONDITIONAL TO: 3-Tesla
The trays have a 2 year expiration date.
MRI Biopsy Tray is **SAFE** in the MR environment when used in the prescribed fashion for all Tesla Strengths.

MODEL # BSM1000-01 **Ship Weight:** 8 lbs
10 TRAYS PER CASE **Ship Dimensions:** 16" x 12" x 12"

NON-FERROMAGNETIC

QTY.	DESCRIPTION
1	5cc syringe luer lock
2	10cc syringe luer lock
1	20cc syringe luer lock
1	20 gauge needle 9cm
1	#11 blade scalpel with plastic handle
1	Plastic forcep fine tip
6	4x4's 12ply
2	2x2's 12ply
1	ChloroPrep® 3ml
1	Betadine swabs (3pk)
1	Sterile pen with ruler
1	Sheet pre-labeled labels
1	Alcohol prep
1	Steri-strip ½"x4" (3 per card)
1	4oz specimen cup
1	Towel, 17" x 20" whit
1	Prep Tray 10"x5"x2"
1	Outer Tray Wrap, 30" x 30"

MRI CERTIFIED #11 DISPOSABLE SCALPEL AND 20G 9CM NEEDLE TRAY KITS

Non-Clinical testing demonstrated that this product is MR conditional and can be used in the MRI environment according to the following conditions:

1. Static magnetic of 3-Tesla or less
2. Highest spatial gradient magnetic field of 720-Gauss/cm or less

Tray Contents:

- 1 CSR Outer Wrap 28 x 28
- 1 Prep Tray 10 x 5 x 2
- 1 MRI "Non-Ferromagnetic" Disposable Titanium #11 Scalpel - This scalpel is safe to use for any interventional procedure for the initial puncture of the skin
- 1 MRI "Non-Ferromagnetic" Needle 20g 9cm - This needle is completely safe for use as a puncture or injection needle.

MR CONDITIONAL TO: 3-Tesla
MODEL # BSM1014-01 1 KIT
MODEL # BSM1015-01 10 KITS
MODEL # BSM1016-01 20 KITS

NON-FERROMAGNETIC



NON MAGNETIC TOOLKITS

Features a gasketed, water-tight and dust-tight, submersible design that is resistant to UV, solvents, corrosion, fungus and impact damage. It also features a molded-in hinge, patent pending trigger release latch system, a comfortable snap-down rubber over-molded cushion grip handle, and an automatic ambient pressure equalization valve.



WATERPROOF CASE




TLM1980-01



TLM1980-01



TLM1981-01



TLM1981-01

MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
TLM1980-01	23 Piece Metric/SAE Toolkit	23 lbs	21" x 17" x 7"
TLM1981-01	25 Piece Metric Toolkit	25 lbs	21" x 17" x 7"

MRI TOOLS



6 PIECE SCREW DRIVER KIT

MODEL # **TLM1510-01**

Standard Screwdriver, 5/16" Tip, 3" Blade, Standard Screwdriver, 7/32" Tip, 6" Blade, Standard Screwdriver, 25/64" Tip, 8" Blade, #1 Phillips Screwdriver, #2 Phillips Screwdriver, #3 Phillips Screwdriver

Ship Weight: 34 lbs

Ship Dimensions: 19" x 12" x 4"



6 PIECE TOOL KIT

MODEL # **TLM1508-01**

Screw Driver - 5/16" Tip, 6 1/2" Blade, Standard Screwdriver, 7/32" Tip, 6" Blade, Standard Screwdriver, 25/64" Tip, 8" Blade, #1 Phillips Screwdriver, #2 Phillips Screwdriver, #3 Phillips Screwdriver, #10 Groove Joint Pliers

Ship Weight: 34 lbs

Ship Dimensions: 19" x 12" x 4"

WATER FIRE EXTINGUISHER

The ultimate extinguisher for class A fires, especially where a potential Class C (electrical) hazard exists. The fine spray from the unique misting nozzle provides safety from electrical shock, greatly enhances the cooling and soaking characteristics of the agent and reduces scattering of the burning materials. Unequivocally the best extinguisher for hospital environments, telecommunications facilities and "clean room" manufacturing facilities.

Includes:

- 2 Bottles of De-Ionized Water
- 1 Fill Adapter



MODEL #	DESCRIPTION	SHIP WEIGHT	SHIP L X W X H
FEM1006-01	1 3/4 Gallon Kit	70 lbs	33" x 22" x 16"
FEM1003-01	2 1/2 Gallon Kit	88 lbs	39" x 22" x 17"

MRI CO2 FIRE EXTINGUISHER

MODEL # **FEM1008-01**

MRI 5 lb Carbon Dioxide(Co2) Fire Extinguisher Model 322NM for Class B Flammable Liquids and is Electrically Non-Conductive.

Wall Hanger Included

All metal Valve Construction

Rust Free Aluminum Cylinder with High Gloss Polyester Powder Paint

5 YEAR WARRANTY

Ship Weight: 14 lbs


Ship Dimensions: 18" x 8" x 5"

THIS EXTINGUISHER IS RATED TO: 7 TESLA




MRI DUST MOP & HANDLE

MODEL # **MTM1146-01** 18" Wide
MODEL # **MTM1147-01** 24" Wide
Ship Weight: 24 lbs
Ship Dimensions: 8" x 8" x 60"



MOP HANDLE

MODEL # **MTM1373-01**
Ship Weight: 26 lbs
Ship Dimensions: 8" x 8" x 65"



MOP BUCKET AND WRINGER

MODEL # **MTM1907-01**
MRI Non-Magnetic Mop Pail with wringer features a divided bucket with a heavy-duty snap-on cone wringer. Wringer slides from side to side accessing both compartments.


Features:
Divided Bucket
Four 2" Casters
15 Quart Bucket

WRINGER SPECIFICATIONS:
Works with up to 16 oz. Mop Heads
Mop Head Available # **MTM1909-01**

Ship Weight: 22 lbs
Ship Dimensions: 14" x 14" x 18"

MRI MOP/DUSTER PADS MR

MODEL # **MTM7012-01**
6" x 20"
MODEL # **MTM7013-01**
6" x 26"
CASE OF 5 EACH
Ship Weight: 2 lbs
Ship Dimensions: 14" x 6" x 4"



WET TUBE MICROFIBER MOP

MODEL # **MTM1918-01**

Features:
Five times the launderability of traditional mops; textured finish loosens and holds dirt until washed
Wide canvas headbands and tailbands are longer lasting and come in a variety of colors for coding
Save money and improve results by converting to microfiber mops
Priced competitively with traditional blended looped yarn mop heads, yet out performs traditional mops in cleaning & durability



PURE HARD SURFACE **DISINFECTANT**

Pure Hard Surface Disinfectant

Pure hard surface disinfectant and food contact sanitizer delivers broad spectrum power to eliminate dangerous germs and foodborne pathogens in food processing, food preparation and public eating places.

FEATURES

- Fast kill times
- GRAS status as contact biocide
- Eliminates bacteria, viruses and fungi
- 24-Hour residual protection
- Odorless, no skin irritation, non-flammable
- No rinse required
- Ideal for use in food processing, food preparation and hospitality

Cleaning Wand & Pads Kit

MTM4550-01 CLEANING KIT INCLUDES:

- 1 wand with disposable pad attached
- 25 additional disposable pads
- 1 32 Oz Pure Hard Surface bottle

MODEL#	DESCRIPTION
MTM4550-01	MRI Cleaning Kit
MTM4551-01	Cleaning Wand with 1 Pad
MTM4552-01	100 Replacement Pads
MTM4553-01	400 Replacement Pads
MTM4554-01	800 Replacement Pads
MTM4508-01	Pure Disinfectant, 32oz Spray Bottle



VELCRO TOP



SPONGE BOTTOM



LADDERS DOUBLE SIDED

375 LB WEIGHT CAPACITY



FIBERGLASS

MODEL #	STYLE	SIZE	SHIP WEIGHT	SHIP L X W X H
Disc.	Aluminum	4 Foot	21 lbs	24" x 9" x 50"
LAM1700-01	Aluminum	5 Foot	28 lbs	26" x 8" x 61"
LAM1702-01	Aluminum	6 Foot	34 lbs	28" x 9" x 74"
LAM1701-01	Aluminum	7 Foot	39 lbs	30" x 8" x 86"
LAM1713-01	Aluminum	8 Foot	41 lbs	32" x 9" x 98"
Disc.	Fiberglass	4 Foot	29 lbs	24" x 9" x 50"
Disc.	Fiberglass	5 Foot	29 lbs	26" x 9" x 62"
Disc.	Fiberglass	6 Foot	37 lbs	28" x 9" x 74"
Disc.	Fiberglass	7 Foot	42 lbs	30" x 9" x 86"
Disc.	Fiberglass	8 Foot	49 lbs	32" x 9" x 98"



ALUMINUM

24" FOLDING SINGLE CLIMB SIDE **STAND**

MODEL #LAM1704-01

24" x 8" Platform, Aluminum

Ship Weight: 34 lbs

Ship Dimensions: 10" x 24" x 27"





MRI BARIATRIC PATIENT **TRANSFER BOARD** MR

MODEL # TRM1011-01

MRI Bariatric Patient Transfer Boards are semi-flexible moving platforms, with hand holes all the way around allowing safe and efficient transfer of patients for general imaging, MRI and CT.

WEIGHT CAPACITY: 300 LBS

MODEL #	DESCRIPTION	SHIP L X W X H
TRM1011-01	Transfer Board	72" x 38" x 5"




MRI ANTI STATIC BARIATRIC **PATIENT SHIFTER** MR

MODEL # TRM1013-01

MRI Anti Static Bariatric Patient Shifters make lateral transfers for patients up to 600 lbs easier and safer and reduces caretaker back strain by minimizing friction on transfer services. Elongated shape brings fourteen handholds to ergonomically correct position. Extra-wide design better supports excess tissue. It even supports patient's head. Flexible plastic is great for uneven surfaces. Clean with hospital-grade disinfectant.

WEIGHT CAPACITY: 250 LBS

MODEL #	DESCRIPTION	SHIP L X W X H
TRM1013-01	Patient Shifter	72" x 38" x 5"



MRI SLIDE BOARD PATIENT **TRANSFER SHEET** MR

MODEL # TRM1014-01

MRI Slide Board Patient Transfer Sheets can be placed under a patient on and off the operating room table in the sitting or supine position.

WEIGHT CAPACITY: 250 LBS

MODEL #	DESCRIPTION	SHIP L X W X H
TRM1014-01	Transfer Sheet	23" x 31" x 5"



MRI TRANSFER DEVICE **STORAGE RACK** MR

MODEL # TRM1021-01

MRI Transfer Device Storage Racks adjust to hold transfer board vertically or horizontally. One board can be stored per rack. Transfer Device Storage Racks are 18" L x 2.5" H x 1" D.

MODEL #	DESCRIPTION	SHIP L X W X H
TRM1021-01	Storage Rack	19" x 5" x 5"

MRI ANTI STATIC CHARCOAL GREY TRANSFER BOARD MR

MODEL # TRM1019-SZ

MRI Anti Static Charcoal Grey Transfer Boards make transfer over uneven surfaces safe and simple. Anti Static Charcoal Grey Transfer Boards easily slide on a bed or stretcher pad. Available either 18" wide or 22" wide. imaging, MRI and CT.

WEIGHT CAPACITY: 250 LBS




MODEL #	DESCRIPTION	SHIP L X W X H
TRM1019-01-SM	Transfer Board	72" x 18" x 5"
TRM1019-01-LG	Transfer Board	72" x 18" x 5"

MRI PLASTIC TRANSFER BOARD MR

MODEL # TRM1016-01

MRI Plastic Transfer Boards bridge the gap between two patient transfer points. The Transfer Board is made of flexible plastic. It is comfortable to sit on and is easily removed after use. Plastic Transfer Boards help the patient transfer between surfaces of similar or slightly different heights.

WEIGHT CAPACITY: 300 LBS




MODEL #	DESCRIPTION	SHIP L X W X H
TRM1016-01	Transfer Board	72" x 22" x 5"

MRI CURVED BANANA TRANSFER BOARD MR

MODEL # TRM1017

MRI Curved Banana Transfer Boards enable a variety of positioning options but allows for the natural sweeping motion of an assisted transfer. The Curved Banana Transfer Board is ideal for bed, car, wheelchair, and toileting transfers. The Curved Banana Transfer Board is made from durable, semi-rigid polyethylene that can be easily wiped down for cleaning and infection control.

WEIGHT CAPACITY: 250 LBS



MODEL #	DESCRIPTION	SHIP L X W X H
TRM1017-01	Transfer Board	30" x 15" x 10"

MRI STATIC TREATED WHITE TRANSFER BOARD MR

MODEL # TRM1020

MRI Static Treated White Transfer Boards make transfers over uneven surfaces safe and simple. Static Treated White Transfer Boards easily slide on a bed or stretcher pad.

WEIGHT CAPACITY: 250 LBS



MODEL #	DESCRIPTION	SHIP L X W X H
TRM-1020-01	Transfer Board	72" x 22" x 5"

2025 PRODUCT GUIDE //

www.clmed.com



SAND BAG POSITIONERS PG.31



MRI STEP STOOL PG.15



NON-FERROMAGNETIC LAUNDRY HAMPER PG.24




NEONATAL NOISE GUARDS PG.33



EYE SHADE SLEEP MASK PG.33



SCAN QR CODE OR
VISIT WWW.CLMED.COM



MR
CONDITIONAL
AT **7T** OR LESS